

## Method for HLA typing

The present invention relates to a method for HLA typing by the unambiguous determination of short DNA sequence elements (2-6 bases) at a given position simultaneously on both parental alleles at a selected number of positions in HLA genes, comprised of the steps for each position of a) hybridising a combination of oligonucleotides (primers) complementary to all known sequence variants to a DNA strand upstream of a given position; b) carrying out a primer extension reaction with at least one of the four dNTP substrates substituted by a terminating analog; c) analysing the products by mass spectrometry, with the resulting masses allowing unambiguous identification of the used primers and the added bases. This method is particularly well suited for DNA-based HLA typing and in combination with a suitable selection of sites tested, it is superior in ease of operation to conventional HLA typing methods.

The most important of the genome projects, the complete sequence of the human genome, is finished. This project reveals the complete sequence of the 3 billion bases and the relative positions of all estimated 30.000 genes in this genome. Having this sequence opens unlimited possibilities for the elucidation of gene function and interaction of different genes. In recent years a systematic effort (SNP consortium) has been underway to identify single nucleotide polymorphisms (SNPs) throughout the human genome and so far several million of these differences between different human beings have been identified (dbSNP contained 5.5 million SNPs in October 2003).

Matrix-assisted laser desorption/ionization time-of-flight mass spectrometry (MALDI) has revolutionized the mass spectrometric analysis of biomolecules (Karas, M. & Hillenkamp, F. *Anal. Chem.* **60**, 2299-2301 (1988)). The field of DNA analysis by mass spectrometry was recently extensively reviewed by Tost and Gut (*Mass Spectrometry Reviews*, **21**, 388-418 (2002)) and Sauer and Gut (*Journal of Chromatography B*, **782**, 73-87, (2002)). MALDI has been applied to the analysis of DNA in variations that range from the analysis of PCR products to approaches using allele-specific termination to single nucleotide primer extension reactions and sequencing (Liu, Y.-H., *et al. Rapid Commun. Mass Spectrom.* **9**, 735-743 (1995);

BEST AVAILABLE COPY

BEST AVAILABLE COPY

Ch'ang, L.-Y., *et al.* *Rapid Commun. Mass Spectrom.* 9, 772-774 (1995); Little, D.P., *et al.* *J. Mol. Med.* 75, 745-750 (1997); Haff, L. & Smirnov, I.P. *Genome Res.* 7, 378-388 (1997), Fei, Z., Ono, T. & Smith, L.M. *Nucleic Acids Res.* 26, 2827-2828 (1998); Ross, P., Hall, L., Smirnov, I. & Haff, L. *Nature Biotech.* 16, 1347-1351 (1998); Ross, P.L., Lee, K. & Belgrader, P. *Anal. Chem.* 69, 4197-4202 (1997); Griffin, T.J., Tang, W. & Smith, L.M. *Nature Biotech.* 15, 1368-1372 (1997); Köster, H., Higgins, G.S & Little, D.P. US Patent 6,043,031). These methods are used to genotype previously identified mutations, SNPs, or insertion/deletions (indels). Spin column purification and/or magnetic bead technology, reversed-phase purification, or ion-exchange resins are frequently applied prior to mass spectrometric analysis.

The GOOD assay (IG Gut et S. Beck: US 6,268,812 ; IG Gut et al: US 6,503,710) is a method for SNP genotyping that uses MALDI mass spectrometry for detection (Sauer et al. 28, e13 and e100 (2000)). Allele-distinction is based on primer extension. In order to make products more amenable to MALDI analysis a substantial part of the primer is removed prior to mass spectrometric analysis. A further element that is included is charge tagging. This means that the final product is conditioned such that it carries either a single positive or a single negative charge. Generally this is achieved by alkylation of a phosphorothioate backbone and in some instances including a quaternary ammonium group to the penultimate base of the primer. The attachment of the quaternary ammonium group gives options for the design of multiplexes - individual SNPs can be moved up or down in the mass spectrum to achieve optimal resolution and separation.

The major histocompatibility complex (MHC) of humans is a cluster of genes on chromosome 6p21. It is of greatest importance as many diseases show association with genes in this region of the genome. All human leukocyte antigen (HLA) coding genes are found in the MHC. The HLA genes are highly variable and implicated in tissue transplantation, immunity and autoimmune disease such as diabetes, psoriasis, lupus, Crohn's disease, colitis, arthritis, and others. The HLA class I genes are HLA-A, HLA-B, HLA-C, ..... The HLA class II genes are HLA-DR, HLA-DQ, HLA-DP,....

HLA typing methods differ dramatically in their approaches. Serological tests can be carried out but have only limited resolution. In the last 15 years the DNA sequence of the MHC has been extensively studied and high resolution typing now makes use of a wealth of DNA sequence information. Methods for DNA based HLA typing range from SSA (sequence specific amplification) where combinations of primers that are specific for different alleles are used to carry out PCR (US 5,545,526). Primers are combined in a way that the sizing of the PCR products allows unambiguous assignment of present base combinations. Multiple combinations are used to identify HLA types. The procedure works its way through a tree of combinations starting with a grouping into rough classes from where on further tests are carried out with specific reagents to subdivide in a class. This method is also known as SSP (sequence specific primers). An alternative method is termed SSOP (sequence specific oligonucleotide probes; US 6,503,707). Here a locus specific PCR is carried out followed by hybridisation with sequence specific oligonucleotide probes. As sequencing technology (and in particular the software for sequence calling) has dramatically improved over the last decade it now is also possible to gain a good degree of identification of HLA types by sequencing (WO 98/35059). Effectively a locus-specific PCR product is sequenced. Problems that arise here are that heterozygous individuals occasionally give rise to ambiguous haplotype calls that can not be resolved (Robinson, J.; Waller, M.J.; Marsh, St.G.E.: "Exon Identities and Ambiguous Typing Combinations"; IMGT/HLA Database; October 2003). The inclusion of allele-specific PCR helps achieve certainty. Resolution requires multiple products per locus to be generated and sequenced. However, as sequencing results can be very convoluted the interpretation in absence of allele-specific PCR can be cumbersome. All together the sequence-based typing requires many iterations in application. Reference strand mediated conformation analysis (RSCA) is a method used to study samples that potentially have a previously unknown sequence in their HLA (Correl et al., Tissue Antigens 56, 82-86, 2000). For a recent review for the reasoning of HLA typing as well as methodological advances see Petersdorf et al. (Tissue Antigens, 61, 1-11, 2003).

The inventors have thus set themselves the task of providing an easy method for the simultaneous capture of all parental mini-haplotypes in highly polymorphic regions of genomes. The procedure has to be executable on a cost-effective genotyping platform. The method should be particularly applicable for HLA typing. It is an aim  
5 to resolve frequent and rare HLA alleles as well as possible.

The object of the present invention is a method for HLA typing by the unambiguous determination of short DNA sequence elements (2-6 bases) simultaneously on both parental alleles at a selected number of positions in HLA genes, comprised of the  
10 steps for each position of a) hybridising a combination of oligonucleotides (primer pool) complementary to all known sequence variants to a DNA strand upstream of a given position; b) carrying out a primer extension reaction with at least one of the four dNTP substrates substituted by a terminating analog; c) analysing the products by mass spectrometry, with the resulting masses allowing unambiguous  
15 identification of the used primers and the added bases.

In the present invention:

- "HLA" means the human leukocyte antigen locus on chromosome 6p21, consisting of HLA genes (HLA-A, HLA-B, HLA-C, HLA-DRB1,...) that are  
20 used to determine the degree of matching, for example, between a recipient and a donor of a tissue graft.
- "HLA typing" means the identification of a known HLA allele of a given locus (HLA-A, HLA-B, HLA-C, HLA-DRB1,...).
- "HLA allele" means a nucleotide sequence within a locus on one of the two  
25 parental chromosomes.
- "HLA-A" means the DNA sequence of exons 2 and 3 of the HLA-A gene.
- "HLA-B" means the DNA sequence of exons 2 and 3 of the HLA-B gene.
- "HLA-DRB1" means the DNA sequence of exon 2 of the HLA-DRB1 gene.
- "Polymorphism" means individual positions in a DNA sequence that exist in  
30 different variants.
- "Haplotype" means the DNA sequence of one of the two alleles in a give region of the genome.

- "Mini-haplotype" means 2-6 contiguous bases on one parental allele.
  - "Primer pools" or "pools of primers" means sets of primers that are used in one primer extension reaction. For each known HLA allele at least one primer is in the pool that is completely complementary in sequence. This assures perfect annealing. Mismatches that are more than 4 bases from the 3' end of the primer do not affect the results of the GOOD assay, as all of those bases are removed by 5'phosphodiesterase after the primer extension reaction. Primers of the pool containing mismatches in the last few bases are not extended by the DNA polymerase and thus not observable.
  - "MALDI mass spectrometer" means a mass spectrometer that uses matrix-assisted laser desorption/ionization for the volatilisation of a sample and time-of-flight analysis for mass separation.
  - "Subgroup" means alleles, which are identical after the mini-haplotyping of the first set of selected positions. For the high resolution typing we resolve subgroups generated with 10 mini-haplotyping reactions. The criteria for resolving subgroups are: a) they still contain alleles with different two-digit types, b) subgroups with more than four alleles, and c) subgroups with frequent alleles (see list below).
- Here we show a methodology for the determination of sequence motifs of 2-6 bases in very polymorphic regions of genomes. In principle this method equates to the determination of mini-haplotypes of 2-6 bases. The individual parental mini-haplotypes can be determined in one reaction without ambiguities. This methodology is applied to a chosen set of positions for HLA typing of HLA-A, HLA-B, and HLA-DRB1. The sets disclosed here have different purposes. First sets of 19, 19, and 10 positions are suggested to distinguish a maximum of HLA alleles in HLA-A, HLA-B, and HLA-DRB1, respectively, with respect to differentiating alleles that are frequent in the general population from ones that are rare. The frequent alleles that were screened for are A\*0101, A\*0201, A\*0301, A\*2301, A\*2402, A\*2902, A\*3001 and A\*3002 for HLA-A, B\*0702, B\*0801, B\*1302, B\*1501, B\*1801, B\*3501, B\*3503, B\*4001, B\*4402, B\*4403, B\*5101 and B\*5701 for HLA-B, and DRB1\*0101, DRB1\*0301, DRB1\*0401, DRB1\*0701,

DRB1\*1101, DRB1\*1104, DRB1\*1302 and DRB1\*1501 for HLA-DRB1. This set of markers provides unambiguous identification of frequent HLA alleles with 93.4 - 100 % certainty in HLA-A, 97.6 - 100 % in HLA-B, and 97.2 - 100 % in HLA-DRB1.

- 5 A second set of 10 positions each in HLA-A, HLA-B, and HLA-DRB1, respectively are described that provide a maximum number of subgroups, that can then be further resolved by the addition of a set of subgroup specific positions. Again the ten positions in each locus were chosen on the basis of providing best distinction between the frequent HLA alleles listed above from the rest of the HLA
- 10 alleles (rare). This resulted in groups containing 2-30 HLA alleles depending on the locus. Within each group a number of positions can be tested to provide resolution between the HLA alleles within the group. The number of positions that have to be additionally analysed range from 1-25 in order to achieve 4-digit resolution. With this technology HLA typing can be carried out at a substantially reduced cost with a
- 15 proven high-throughput detection platform (MALDI mass spectrometry).

In a preferred embodiment of the method of the invention, the DNA strand of step a) is produced by a DNA replication procedure such as PCR or rolling circle replication.

- 20 A set of locus-specific PCR reactions for the selective amplification of each locus is described by the International Histocompatibility Working Group, Technical Manuals ([www.ihwg.org/tmanual/Tmcontents.htm](http://www.ihwg.org/tmanual/Tmcontents.htm)).

In a very preferred embodiment of the method of the invention, a combination of primers (pools of primers) contains slightly varying sequences so that all known

25 sequences of the HLA alleles are accommodated by a perfectly matching primer.

The pool of primers guarantees that at least one primer is perfectly matched. The hybridised oligonucleotides of the primer pool are extended onto a polymorphic position. A requirement is that the added base together with the base composition of the primer gives a unique mass. The detection of this mass in the mass

30 spectrometric profile indicates the presence of a sequence containing both the complementary sequence of the primer and the added base. In order to make all primers of a primer pool distinguishable by mass it is possible to add different mass

shifting agents to the primers. The easiest way to accomplish this is by using charge/mass tagging technology such as is used in the GOOD assay. The penultimate base from the 3' end of the primer is amino-modified and used to add tags via NHS-ester chemistry. The pools of primers of course contain primers that sometimes differ by as little as one base. Sequences identical in base content can still be distinguished by the suitable selection of mass tags. Also, we have found that a primer carrying a mismatch in the last eight bases from the 3' end even if it anneals is not extended by the polymerase and thus screened out. This might be due to insufficient hybridisation or a resistance of the DNA polymerase to attach or extend when a mismatch is present. We thus make use of two effects for our mini-haplotyping: 1) allele-specific hybridisation and 2) allele-specific primer extension. Mismatches that are further than four bases away from the 3' end of the extension primer do not result in increased complexity of the mass spectra as they are removed in the 5'phosphodiesterase digestion step of the GOOD assay.

In a preferred embodiment of the method of the invention, mass shifting tags are added to the individual primers sequences of a primer pool to make them uniquely distinguishable once the terminating base is added.

In another preferred embodiment of the method of the invention, termination products for known alleles are generated by extending the perfectly hybridised primer with a combination of dNTPs and ddNTPs or analogues thereof with a DNA polymerase to generate specific termination products to make them uniquely distinguishable by their mass.

In a preferred embodiment of the method of the invention, the GOOD assay is used. It typically applies single base primer extension, thus only the four terminating bases (ddNTPs) or synthetic analogues with the same qualities in terms of DNA polymerase tolerance are used for primer extension.  $\alpha$ -S-ddNTPs are very suitable analogues.

In a preferred embodiment of the method of the invention, mass spectrometry, in particular MALDI or ESI mass spectrometry is used for analysis of the masses of products.

For HLA typing a set of said mini-haplotyping assays has to be carried out to achieve sufficient information content.

For HLA typing of HLA-A the preferred set of assays are those of positions 98, 414, 539, 282, 571, 368, 256, 292, 238, 270, 453, 527, 502, 81, 268, 559, 92, 123 and 396 (according to the numbering of the HLA-A gene starting at cDNA sequence position 1 of exon 1; see Figure 1). This results in medium resolution HLA typing. The input criteria for the selection are the frequency of HLA alleles.

Some HLA types are identified unambiguously.

For HLA typing of HLA-B accordingly the following positions are preferably analysed by mini-haplotyping assays to achieve medium resolution: 539, 419, 559, 412, 272, 362, 302, 363, 206, 369, 259, 97, 583, 292, 222, 527, 418, 435 and 571 (according to the numbering of the HLA-B gene starting at cDNA sequence position 1 of exon 1; see Figure 2).

For HLA typing of HLA-DRB1 accordingly the following positions are preferably analysed by mini-haplotyping to achieve medium resolution: 125, 196, 197, 227, 261, 286, 299, 308, 341 and 345 (according to the numbering of the HLA-DRB1 gene starting at cDNA sequence position 1 of exon 1; see Figure 3).

In a preferred embodiment for high resolution HLA typing of HLA-A positions 98, 414, 539, 282, 571, 368, 256, 292, 238 and 270 (according to the numbering of the HLA-A gene starting at cDNA sequence position 1 of exon 1; see Figure 4) are used for mini-haplotyping to generate sub-groups (HLA-A\_A, HLA-A\_B, HLA-A\_C, HLA-A\_D, HLA-A\_E, HLA-A\_F, HLA-A\_G, HLA-A\_H, HLA-A\_I, HLA-A\_J, HLA-A\_K, HLA-A\_L, HLA-A\_M, HLA-A\_N, and HLA-A\_O; see Table I).

Positions 224, 268, 376, 502, 561 and 616 are preferably analysed to resolve subgroup HLA-A\_A (sequences identical over exons 2 and 3 for alleles A\*29010101 and A\*29010102); positions 126 and 526 to resolve subgroup HLA-A\_B; positions 81, 90, 92, 212, 214, 257, 265, 299, 302, 404, 420, 427, 453, 485, 489 and 502 to resolve subgroup HLA-A\_C (sequences identical over exons 2 and 3 for alleles A\*24020101, A\*24020102L, A\*240203, A\*2409N and A\*2411N); positions 160, 200, 362 and 524 to resolve subgroup HLA-A\_D; positions 180, 299, 301, 302, 346, 418, 453, 517, 524, 526, 527, 557, 559 and 560 to resolve subgroup HLA-A\_E; positions 299, 301, 302, 341 and 583 to resolve subgroup HLA-A\_F; positions 127, 341, 399, 480, 502, 503, 524, 526, 527, 553, 559, 560 and 565 to resolve subgroup HLA-A\_G; positions 228, 233, 463, 519, 530 and 583 to resolve



subgroup HLA-A\_H; positions 102, 275, 317, 362, 418, 419, 497, 524, 555, 595 and 618 to resolve subgroup HLA-A\_I (sequences identical over exons 2 and 3 for alleles A\*680102 and A\*6811N); positions 92, 331, 453, 524, 559, 560 and 564 to resolve subgroup HLA-A\_J; positions 78, 81, 123, 125, 142, 144, 194, 268, 294, 324, 355, 362, 396, 403, 419, 453, 456, 477, 493, 517, 524, 526, 527, 559 and 560 to resolve subgroup HLA-A\_K (sequences identical over exons 2 and 3 for alleles A\*02010101, A\*02010102, A\*020108, A\*0209, A\*0243N and A\*0266); positions 113, 299, 301, 302, 308, 311, 523, 524 to resolve subgroup HLA-A\_L; positions 171, 363, 498 and 559 to resolve subgroup HLA-A\_M; positions 376, 426, 527, 555, 557 and 595 to resolve subgroup HLA-A\_N; position 299 to resolve subgroup HLA-A\_O.

**TABLE I**

Subgroups of HLA-A	Alleles of Subgroups	Positions to resolve Subgroups
HLA-A_A	A*29010101, A*29010102, A*290201, A*290202, A*2904, A*2906, A*2908N, A*2909	224, 268, 376, 502, 561, 616
HLA-A_B	A*3002, A*3009, A*3012	126, 526
HLA-A_C	A*24020101, A*24020102L, A*240202, A*240203, A*240204, A*2404, A*2405, A*2408, A*2409N, A*2411N, A*2420, A*2421, A*2425, A*2426, A*2427, A*2429, A*2432, A*2435, A*2436N, A*2437, A*2438, A*2439	81, 90, 92, 212, 214, 257, 265, 299, 302, 404, 420, 427, 453, 485, 485, 489, 502
HLA-A_D	A*0206, A*0214, A*0221, A*0251, A*0257	160, 200, 362, 524
HLA-A_E	A*250101, A*250102, A*2601, A*2604, A*2605, A*2609, A*2610, A*2611N, A*2612, A*2614, A*2615, A*2617, A*2618, A*6603	180, 299, 301, 302, 346, 418, 453, 517, 524, 526, 527, 557, 559, 560
HLA-A_F	A*2502, A*2613, A*6601, A*6602, A*6604	299, 301, 302, 341, 583
HLA-A_G	A*110101, A*110102, A*1102, A*1103, A*1104, A*1105, A*1107, A*1109, A*1112, A*1113, A*1114, A*1115	127, 341, 399, 480, 502, 503, 524, 526, 527, 553, 559, 560, 565
HLA-A_H	A*3301, A*330301, A*330302, A*3304, A*3305, A*3306, A*3307	228, 233, 463, 519, 530, 583
HLA-A_I	A*680101, A*680102, A*680103, A*6807, A*6811N, A*6812, A*6816, A*6817, A*6819, A*6821, A*6822, A*6823, A*6824	102, 275, 317, 362, 418, 419, 497, 524, 555, 595, 618
HLA-A_J	A*2301, A*2303, A*2305, A*2306, A*2307N, A*2308N, A*2310, A*2413	92, 331, 453, 524, 556, 560, 564
HLA-A_K	A*02010101, A*02010102, A*020102, A*020103, A*020104, A*020105, A*020106, A*020107, A*020108, A*020109, A*0204, A*0209, A*0216, A*0224, A*0225, A*0226, A*0229, A*0230, A*0231, A*0232N, A*0240, A*0242, A*0243N, A*0258, A*0259, A*0260, A*0264, A*0266, A*0267, A*0253N	78, 81, 123, 125, 142, 144, 194, 268, 294, 324, 355, 362, 396, 403, 419, 453, 419, 453, 456, 477, 493, 517, 524, 526, 527, 559, 560
HLA-A_L	A*3201, A*3203, A*3206, A*7401, A*7402, A*7403, A*7408, A*7409	113, 299, 301, 302, 308, 311, 523, 524
HLA-A_M	A*010101, A*010102, A*0103, A*0104N, A*0108, A*0109	171, 363, 498, 559
HLA-A_N	A*03010101, A*03010102, A*0303N, A*0304, A*0305, A*0306, A*0307, A*0311N	376, 426, 527, 555, 557, 595
HLA-A_O	A*2504, A*2608	299

In a preferred embodiment for high resolution, HLA typing of HLA-B positions 539, 419, 559, 412, 272, 362, 302, 363, 206 and 369 (according to the numbering of the HLA-B gene starting at cDNA sequence position 1 of exon 1; see Figure 5) are used for mini-haplotyping to generate sub-groups (HLA-B\_A, HLA-B\_B, HLA-B\_C, HLA-B\_D, HLA-B\_E, HLA-B\_F, HLA-B\_G, HLA-B\_H, HLA-B\_I, HLA-B\_J, HLA-B\_K, HLA-B\_L, HLA-B\_M, HLA-B\_N, HLA-B\_O, HLA-B\_P, HLA-B\_Q, HLA-B\_R, HLA-B\_S, HLA-B\_T, HLA-B\_U, HLA-B\_V, HLA-B\_W, HLA-B\_X, HLA-B\_Y, HLA-B\_Z, HLA-B\_AA, HLA-B\_AB and HLA-B\_AC; see Table II).

Positions 259, 341 and 473 are preferably analyzed to resolve subgroup HLA-B\_A (sequences identical over exons 2 and 3 for alleles B\*0801 and B\*0819N); positions 106, 144, 222, 259, 273, 311, 313, 418, 445, 493, 528 and 540 to resolve subgroup HLA-B\_B (sequences identical over exons 2 and 3 for alleles B\*44020101, B\*44020102, B\*4419N and B\*4427); positions 319, 416, 545 and 572 to resolve subgroup HLA-B\_C; positions 106, 131, 165, 215, 243, 277, 292, 322, 481, 582, 603 and 616 to resolve subgroup HLA-B\_D; positions 106, 146, 165, 181, 238, 259, 263, 292, 328.1/329(insert for B\*1579N), 379, 435, 453, 463, 485, 526, 571, 572 and 583 to resolve subgroup HLA-B\_E (sequences identical over exons 2 and 3 for alleles B\*15010101 and B\*15010102); positions 142, 171, 255, 257, 395, 430, 544, 566 and 572 to resolve subgroup HLA-B\_F; positions 117, 247, 248, 277, 345, 418, 489 and 527 to resolve subgroup HLA-B\_G (sequences identical over exons 2 and 3 for alleles B\*270502, B\*270504 and B\*2713); positions 134, 141, 200, 213, 259, 304 and 527 to resolve subgroup HLA-B\_H; positions 83, 141, 211, 222, 242, 322, 404, 414, 435, 463, 502, 527, 544, 571, 572 and 583 to resolve subgroup HLA-B\_I (sequences identical over exons 2 and for alleles B\*510101, B\*510105, B\*5111N, B\*5130 and B\*5132); positions 103, 142, 222, 243, 259, 292, 477, 486 and 499 to resolve subgroup HLA-B\_J (sequences identical over exons 2 and 3 for alleles B\*400101 and B\*400102); positions 103, 259, 292, 295, 527 and 583 to resolve subgroup HLA-B\_K (sequences identical over exons 2 and 3 for alleles B\*180101 and B\*1817N); positions 320 and 500 to resolve subgroup HLA-B\_L; positions 311, 527 and 583 to resolve subgroup HLA-B\_M; positions 119, 292, 259, 319, 425, 527, 546 and 583 to resolve subgroup HLA-B\_N (sequences identical over exons 2 and 3 for alleles B\*350101, B\*3540N and B\*3542); positions 97, 142, 245 and 527 to resolve subgroup HLA-B\_O; positions 97 and 175 to resolve subgroup HLA-B\_P; positions

**TABLE II**

<i>Subgroups of</i>	<i>Alleles of the subgroup</i>	<i>Positions to resolve</i>
<i>HLA-B</i>		<i>Subgroups</i>
HLA-B_A	B*0801, B*0808N, B*0810, B*0818, B*0819N	259, 341, 473
HLA-B_B	B*44020101, B*44020102S, B*440202, B*440203, B*4405, B*4411, B*4412, B*4419N, B*4422, B*4423N, B*4424, B*4425, B*4427, B*4433, B*4434, B*4435	106, 144, 222, 259, 273, 311, 313, 418, 445, 493, 528, 540
HLA-B_C	B*4415, B*4501, B*4503, B*4504, B*4505	319, 416, 545, 572
HLA-B_D	B*070201, B*070202, B*070203, B*070204, B*0703, B*0716, B*0721, B*0722, B*0723, B*0729, B*0730, B*0733, B*0735	106, 131, 165, 215, 243, 277, 292, 322, 481, 582, 603, 616
HLA-B_E	B*15010101, B*15010102, B*150102, B*150103, B*150104, B*1512, B*1514, B*1515, B*1519, B*1528, B*1533, B*1534, B*1538, B*1560, B*1570, B*1571, B*1575, B*1578, B*1579N, B*1581, B*1582	106, 146, 165, 181, 238, 259, 263, 292, 328.1/329, 379, 435, 453, 463, 485, 526, 571, 572, 583
HLA-B_F	B*440301, B*4413, B*4426, B*4429, B*4430, B*4432, B*4436, B*4437, B*4438, B*4439	142, 171, 255, 257, 395, 430, 544, 566, 572
HLA-B_G	B*2703, B*270502, B*270503, B*270504, B*270505, B*270506, B*2709, B*2710, B*2713, B*2716, B*2717	117, 247, 248, 277, 345, 418, 489, 527
HLA-B_H	B*5107, B*520101, B*520102, B*520103, B*520104, B*5203, B*5204, B*5205	134, 141, 200, 213, 259, 304, 527
HLA-B_I	B*510101, B*510102, B*510103, B*510104, B*510105, B*510201, B*510202, B*5103, B*5109, B*5111N, B*5112, B*5114, B*5118, B*5119, B*5123, B*5124, B*5126, B*5127N, B*5128, B*5130, B*5132, B*5133	83, 141, 211, 222, 242, 322, 404, 414, 435, 463, 502, 527, 544, 571, 572, 583
HLA-B_J	B*400101, B*400102, B*400103, B*4010, B*4011, B*401401, B*401402, B*401403, B*4022N, B*4025, B*4043	103, 142, 222, 243, 259, 292, 477, 486, 499
HLA-B_K	B*180101, B*180102, B*1803, B*1804, B*1805, B*1811, B*1812, B*1815, B*1817N	103, 259, 292, 295, 527, 583
HLA-B_L	B*570101, B*5706, B*5708	320, 500
HLA-B_M	B*3527, B*5301, B*5302, B*5306, B*5308	311, 527, 583
HLA-B_N	B*350101, B*350102, B*3507, B*3510, B*3511, B*3521, B*3524, B*3529, B*3540N, B*3541, B*3542, B*5305	119, 292, 259, 319, 425, 527, 546, 583
HLA-B_O	B*5501, B*5502, B*5505, B*5510, B*5516	97, 142, 245, 527
HLA-B_P	B*5401, B*5402, B*5507	97, 175

HLA-B_Q	B*3910, B*670101, B*670102	246, 277
HLA-B_R	B*3803, B*390201, B*390202, B*3913, B*3923	246, 292, 311, 503
HLA-B_S	B*3801, B*380201, B*380202, B*3804, B*3805, B*3809	103, 261, 309, 311, 474
HLA-B_T	B*390101, B*390103, B*390104, B*3904, B*3905, B*3912, B*3922, B*3925N, B*3926	97, 103, 106, 243, 259, 292, 404, 524
HLA-B_U	B*3503, B*3513, B*3536	259, 320
HLA-B_V	B*0734, B*5504	106
HLA-B_W	B*4047, B*4431	97
HLA-B_X	B*4002, B*4027, B*4029, B*4035, B*4040, B*4045	97, 106, 257, 418, 463
HLA-B_Y	B*400104, B*4004	106
HLA-B_Z	B*4012, B*4046, B*4803	106, 144
HLA-B_AA	B*2703, B*270502, B*270503, B*270504, B*270505, B*270506, B*2709, B*2710, B*2713, B*2716, B*2717	117, 247, 248, 283, 345, 418, 489, 527
HLA-B_AB	B*1562, B*4802	106
HLA-B_AC	B*1302, B*1308	548

246 and 277 to resolve subgroup HLA-B\_Q; positions 246, 292, 311 and 503 to resolve subgroup HLA-B\_R; positions 103, 261, 309, 311 and 474 to resolve subgroup HLA-B\_S; positions 97, 103, 106, 243, 259, 292, 404 and 524 to resolve subgroup HLA-B\_T (sequences identical over exons 2 and 3 for alleles B\*390101 and B\*390103); positions 259 and 320 to resolve subgroup HLA-B\_U; position 106 to resolve HLA-B\_V; positions 97 to resolve HLA-B\_W; positions 97, 106, 257, 418 and 463 to resolve HLA-B\_X; position 106 to resolve HLA-B\_Y; positions 106 and 144 to resolve HLA-B\_Z; positions 117, 247, 248, 283, 345, 418, 489, and 527 to resolve HLA-B\_AA; positions 106 to resolve HLA-B\_AB; positions 548 to resolve HLA-B\_AC.

In a preferred embodiment, the method for HLA typing resolves groups A-P of HLA-DRB1.

For high resolution, HLA typing of HLA-DRB1 positions are: 125, 196, 197, 227, 261, 286, 299, 308, 341 and 345 (according to the numbering of the HLA-DRB1 gene starting at DNA sequence position 1 of exon 1; see Figure 6) are used for mini-haplotyping to generate sub-groups (HLA-DRB1\_A, HLA-DRB1\_B, HLA-DRB1\_C, HLA-DRB1\_D, HLA-DRB1\_E, HLA-DRB1\_F, HLA-DRB1\_G, HLA-DRB1\_H, HLA-DRB1\_I, HLA-DRB1\_J, HLA-DRB1\_K, HLA-DRB1\_L, HLA-DRB1\_M, HLA-DRB1\_N, HLA-DRB1\_O, HLA-DRB1\_P; see Table III).

In a very preferred embodiment, positions 123, 174, 250, 278 and 317 are analysed to resolve subgroup HLA-DRB1\_A; positions 192, 203, 256 and 259 to resolve subgroup HLA-DRB1\_B; 256, 260, 317 and 351 to resolve subgroup HLA-DRB1\_C; positions 155, 204, 233, 239, 256, 304, 357 and 366 to resolve subgroup HLA-DRB1\_D; positions 122, 171, 257 and 317 to resolve subgroup HLA-DRB1\_E; positions 164, 167, 171, 230, 235, 306, 317, 321 and 337 to resolve subgroup HLA-DRB1\_F; positions 164, 257, 266 and 303 to resolve subgroup HLA-DRB1\_G; positions 164, 181, 188, 220, 229, 256, 266, 317 and 318 to resolve subgroup HLA-DRB1\_H; position 257 to resolve subgroup HLA-DRB1\_I; positions 181, 239 and 357 to resolve subgroup HLA-DRB1\_J; positions 122, 144, 239, 303, 317, 318 and 321 to resolve subgroup HLA-DRB1\_K (sequences identical over exons 2 and 3 for alleles DRB1\*110101 and DRB1\*110102); positions 118, 161, 257, 260, 318 and 321 to resolve subgroup HLA-DRB1\_L; positions 165, 257, 293 and 303 to resolve subgroup HLA-DRB1\_M (sequences identical over exons 2 and 3 for alleles DRB1\*120101 and DRB1\*1206); positions 177, 240, 256, 257 and 357 to resolve subgroup HLA-DRB1\_N; positions 150 175, 230, 236 and 321 to resolve subgroup HLA-DRB1\_O (sequences identical over exons 2 and 3 for alleles DRB1\*150101 and DRB1\*1513); positions 115, 220 and 317 to resolve subgroup HLA-DRB1\_P.

Another object of the invention is a kit to carry out the procedure. It consists of pooled combinations of primers. The primers that are used in the pools for HLA-A, HLA-B, and HLA-DRB1 and the masses of the genotyping products are listed in Tables IV, V, and VI respectively. CT refers to the mass shifting mass tag that is attached to that primer of the pool.

Another object of the invention is the use of the method of the invention for screening of tissue donors.

In a preferred embodiment, the use is for bone marrow donors in registries for screening of frequent and rare HLA types.

Still another object of the invention is the use of the primers represented in Table IV, V and VI to carry out HLA typing.

**TABLE III**

Subgroups of HLA-DRB1	Alleles of Subgroups	Positions to resolve Subgroups
HLA-DRB1_A	DRB1*070101, DRB1*070102, DRB1*0703, DRB1*0704, DRB1*0705, DRB1*0707	123, 174, 250, 317
HLA-DRB1_B	DRB1*040101, DRB1*040102, DRB1*0409, DRB1*0426, DRB1*0433	192, 203, 256, 259
HLA-DRB1_C	DRB1*0404, DRB1*0410, DRB1*0423, DRB1*0440, DRB1*0444	256, 260, 317, 351.
HLA-DRB1_D	DRB1*040501, DRB1*040502, DRB1*040503, DRB1*040504, DRB1*0408, DRB1*0429, DRB1*0430, DRB1*0445, DRB1*0448	155, 204, 233, 239, 256, 304, 357, 366
HLA-DRB1_E	DRB1*1402, DRB1*1409, DRB1*1413, DRB1*1446, DRB1*1447, DRB1*1448	122, 171, 257, 317
HLA-DRB1_F	DRB1*130101, DRB1*130102, DRB1*130103, DRB1*1315, DRB1*1327,	164, 167, 171, 230, 235, 306, 317, 321, 337
HLA-DRB1_G	DRB1*130201, DRB1*130202, DRB1*1331, DRB1*1339, DRB1*1341	164, 257, 266, 303
HLA-DRB1_H	DRB1*030101, DRB1*030102, DRB1*0307, DRB1*0312, DRB1*0313, DRB1*0315, DRB1*0316, DRB1*0318, DRB1*0322, DRB1*0323	164, 181, 188, 220, 229, 256, 266, 317, 318
HLA-DRB1_I	DRB1*1137, DRB1*1425	257
HLA-DRB1_J	DRB1*110401, DRB1*110402, DRB1*1143, DRB1*1146	181, 239, 357
HLA-DRB1_K	DRB1*110101, DRB1*110102, DRB1*110103, DRB1*110104, DRB1*110105, DRB1*112701, DRB1*112702, DRB1*1130, DRB1*1139	122, 144, 239, 303, 317, 318, 321
HLA-DRB1_L	DRB1*1117, DRB1*140101, DRB1*140102, DRB1*1408, DRB1*1426, DRB1*1438, DRB1*1439	118, 161, 257, 260, 318, 321
HLA-DRB1_M	DRB1*120101, DRB1*120102, DRB1*1206, DRB1*1207, DRB1*1208, DRB1*1209	165, 257, 293, 303
HLA-DRB1_N	DRB1*080101, DRB1*080102, DRB1*080201, DRB1*080202, DRB1*080203, DRB1*0807, DRB1*0811	177, 240, 256, 257, 357
HLA-DRB1_O	DRB1*150101, DRB1*150103, DRB1*150105, DRB1*1503, DRB1*1506, DRB1*1509, DRB1*1513	150 175, 230, 236, 321
HLA-DRB1_P	DRB1*010101, DRB1*0105, DRB1*0107, DRB1*0111	115, 220, 317

16  
TABLE IV

No.	Name	Sequence	CT	Primer Masses	A	C	G	T
1	HLAA 811 1f20	TGCTCGCCCCCAGGCTCCspC <sup>A</sup> spA	0	1098,1	1425,1	1401,3	-	-
2	HLAA 812 1f20	TGCTCGCCCCCAGGCTCTspC <sup>A</sup> spA	0	1113,1	-	1416,3	1452,4	-
3	HLAA 921 1f20	AGGCTCCCACTCCATGAGspC <sup>A</sup> spT	0	1129,1	1456,4	-	-	-
4	HLAA 922 1f20	AGGCTCCCAMTCCATGAGspG <sup>A</sup> spT	0	1169,1	1496,4	-	1512,4	-
5	HLAA 923 1f20	AGGCTCTCASTCCATGAGspG <sup>A</sup> spT	0	1169,1	1496,4	-	1512,4	-
6	HLAA 981 1f20	CCACTCCATGAGGTATTTspC <sup>A</sup> spA	0	1113,1	-	1416,3	-	-
7	HLAA 982 1f20	CCACTCCATGAGGTATTTspC <sup>A</sup> spT	0	1104,1	1431,4	1407,3	-	1422,3
8	HLAA 1231 2r20	GCGATGAAGCGGGGCTCspCspT <sup>A</sup> spC	0	1510,5	-	-	1853,8	-
9	HLAA 1232 2r20	GCGATGAAGCGGGGCTCspTspC <sup>A</sup> spC	-28	1380,4	1707,7	-	-	-
10	HLAA 1233 2r20	GCGATGAAGCGGGGCTTspCspC <sup>A</sup> spC	0	1408,4	-	-	1751,6	-
11	HLAA 1234 2r20	GMGATGAAGCGGGGCTCspCspC <sup>A</sup> spC	0	1393,4	1720,7	-	1736,7	-
12	HLAA 2381 2r20	CTSGTCCCAATACTCCGspGspA <sup>A</sup> spC	0	1497,4	-	1800,6	-	-
13	HLAA 2382 2r20	CYCGTCCCAATACTCCGspGspA <sup>A</sup> spC	0	1497,4	-	1800,6	-	-
14	HLAA 2383 2r20	CTCGTCCCAATACTCCGspGspC <sup>A</sup> spT	0	1488,4	-	1791,6	-	1806,4
15	HLAA 2384 2r20	CTSGTCCCAATACTCAGspGspC <sup>A</sup> spC	0	1473,4	-	1776,6	-	-
16	HLAA 2385 2r20	CYGGTCCCAATACTCCGspGspC <sup>A</sup> spC	0	1473,4	-	1776,6	-	-
17	HLAA 2386 2r20	CMGGTCCCAATACTCCGspGspC <sup>A</sup> spC	0	1473,4	-	1776,6	-	-
18	HLAA 2387 2r20	CYCGTCCCAATACTCCGspGspC <sup>A</sup> spC	0	1473,4	-	1776,6	-	-
19	HLAA 2561 1r19	CTTCATATTCGGTGTCTCspC <sup>A</sup> spT	0	1089,1	-	1392,3	1432,4	-
20	HLAA 2562 1r19	CTTCACWITCCGTGTCTCspC <sup>A</sup> spT	0	1089,1	-	1392,3	1432,4	-
21	HLAA 2563 1r19	CTTCACATKCCGTGTCTGspC <sup>A</sup> spA	0	1138,1	-	-	1481,4	-
22	HLAA 2564 1r19	CTTCACATTCGGTGTGTspC <sup>A</sup> spC	0	1089,1	-	-	1432,1	-
23	HLAA 2565 1r19	CYTACATTCGGTGTGTspC <sup>A</sup> spC	0	1089,1	-	-	1432,1	-
24	HLAA 2566 1r19	CTTCACRTTCGGTGTCTCspC <sup>A</sup> spC	0	1074,1	-	1377,3	1417,4	-
25	HLAA 2567 1r19	CTTCASTTGCCGTGTCTCspC <sup>A</sup> spC	0	1074,1	-	1377,3	1417,4	-
26	HLAA 2568 1r19	CTTCAGTTKCCGTGTCTCspC <sup>A</sup> spC	0	1074,1	-	1377,3	1417,4	-
28	HLAA 2681 1f20	ATTGGGACCGGAACACACspG <sup>A</sup> spG	0	1154,1	1481,4	1457,3	-	-
29	HLAA 2682 1f20	ATTGGGACCTGCAGACACspG <sup>A</sup> spG	0	1154,1	1481,4	1457,3	-	-
30	HLAA 2683 1f20	ATTGGGACsAGGAGACACspG <sup>A</sup> spG	0	1154,1	1481,4	1457,3	-	-
31	HLAA 2684 1f20	ATTGGGACSGGGAGACACspG <sup>A</sup> spG	0	1154,1	1481,4	1457,3	-	-
32	HLAA 2685 1f20	ATTGGGACsAGGAGACAGspG <sup>A</sup> spG	0	1194,1	1521,4	-	-	-
33	HLAA 2701 1r19	CTGTGAGTGGGCCTTCspA <sup>A</sup> spT	0	1113,1	1440,4	-	-	-
34	HLAA 2702 1r19	CTGTGACTGGGCCYTspA <sup>A</sup> spC	-14	1084,1	1411,4	-	1427,4	1402,4
35	HLAA 2703 1r19	CTGTGAGTGGSCCTTCspA <sup>A</sup> spC	-14	1084,1	1411,4	-	1427,4	1402,4
36	HLAA 2821 1f20	ACACGGAATGTGARGGGCspC <sup>A</sup> spA	0	1098,1	-	1401,3	1441,3	-
37	HLAA 2822 1f20	ACASGGAAGTGAAGGCCspC <sup>A</sup> spA	0	1098,1	-	1401,3	1441,3	-
38	HLAA 2823 1f20	ACACGGCAWGTGAAGGCCspC <sup>A</sup> spA	0	1098,1	-	1401,3	1441,3	-
39	HLAA 2824 1f20	ACACGGAACGTGAAGGCCspC <sup>A</sup> spA	0	1098,1	-	1401,3	1441,3	-
40	HLAA 2825 1f20	ACACGGAATRTGAAGGCCspC <sup>A</sup> spA	0	1098,1	-	1401,3	1441,3	-
41	HLAA 2921 2f20	TGAAGGCCCACTCACAGspAspG <sup>A</sup> spT	-14	1498,4	-	1801,6	-	-
42	HLAA 2922 2f20	TGAAGGCCCACTCACAGspGspC <sup>A</sup> spT	0	1488,4	-	-	1831,7	-
43	HLAA 2923 2f20	TGAAGGCCCACTCACAGspAspT <sup>A</sup> spT	0	1589,6	-	-	1932,9	-
44	HLAA 2924 2f20	TGARGGCCAGTCACAGspAspC <sup>A</sup> spT	0	1427,4	-	1775,6	1815,7	-
45	HLAA 2925 2f20	TGAAGGCCCACTCACAGspAspC <sup>A</sup> spT	0	1427,4	-	1775,6	1815,7	-
46	HLAA 3681 1f20	TCACACCATCCAGATAATspG <sup>A</sup> spC	0	1129,1	1456,4	-	-	-
47	HLAA 3682 1f20	TCACACCATCCAGMTAATspG <sup>A</sup> spT	0	1144,1	1471,6	1447,1	1487,4	1462,3
48	HLAA 3683 1f20	TCACACCSTCCAGAGGATspG <sup>A</sup> spT	0	1144,1	1471,6	1447,1	1487,4	1462,3
49	HLAA 3684 1f20	TCACACCVTCCAGATGATspG <sup>A</sup> spT	0	1144,1	1471,6	1447,1	1487,4	1462,3
50	HLAA 3961 2r20	GCTGGTACCCGCGGAGspGspA <sup>A</sup> spG	0	1537,4	-	-	1880,7	-



51	HLAA 3962 2r20	GCCGGTACCCGCGGAGspTspA <sup>spA</sup>	0	1496,4	-	-	1839,7	-
52	HLAA 3963 2r20	GGTGGTACCCGYGCAGspGspA <sup>spA</sup>	0	1496,4	-	-	1839,7	-
53	HLAA 3964 2r20	GGTGGTACCCGCGAGspGspA <sup>spA</sup>	0	1521,5	-	-	1864,8	1839,7
54	HLAA 3965 2r20	GTTCATACCCGCGGAGspGspA <sup>spA</sup>	0	1521,5	-	-	1864,8	1839,7
55	HLAA 3966 2r20	GSTGGTACCCGCGGAGspGspA <sup>spA</sup>	0	1521,5	-	-	1864,8	1839,7
56	HLAA 3967 2r20	GCCGGTACCCGCGGAGspGspA <sup>spA</sup>	0	1521,5	-	-	1864,8	1839,7
57	HLAA 4141 1f20	CGCTTCCTCCGCGGGTATspG <sup>spA</sup>	0	1153,1	1480,1	-	-	-
58	HLAA 4142 1f20	CGCTTCCTCTGCGGGTACspC <sup>spA</sup>	0	1098,1	-	1401,3	1441,4	-
59	HLAA 4143 1f20	CGCTTCCTGCGCGGGTACspC <sup>spA</sup>	0	1098,1	-	1401,3	1441,4	-
60	HLAA 4144 1f20	CGCTTCCTCCACGGGTACspC <sup>spA</sup>	0	1098,1	-	1401,3	1441,4	-
61	HLAA 4145 1f20	CGMTTCCTCCGCGGGTACspC <sup>spA</sup>	0	1098,1	-	1401,3	1441,4	-
62	HLAA 4146 1f20	CGCCTCCTCCGCGGGTACspC <sup>spA</sup>	0	1098,1	-	1401,3	1441,4	-
63	HLAA 4147 1f20	CACTTCCTCCGCGGGTACspC <sup>spG</sup>	0	1114,1	-	-	1457,4	-
64	HLAA 4148 1f20	CGCTTMCTCCGCGGGTACspC <sup>spG</sup>	0	1114,1	-	-	1457,4	-
65	HLAA 4531 1r20	GTCCAAGAGCGCAGGTCTspT <sup>spC</sup>	0	1206,2	-	-	-	1524,4
66	HLAA 4532 1r20	GTCCAAGAGCGCAGGTCCspT <sup>spC</sup>	0	1191,2	-	-	1534,5	1509,4
67	HLAA 4533 1r20	GTCCAGGAGCTCAGGTCCspT <sup>spC</sup>	0	1191,2	-	-	1534,5	1509,4
68	HLAA 5021 2r20	GGCCGYCTCCCACTTGTspGspC <sup>spT</sup>	0	1463,4	-	-	-	1781,6
69	HLAA 5022 2r20	GGCYGCCTCCCACTTGCspGspC <sup>spT</sup>	0	1448,4	-	1751,6	1791,7	1766,6
70	HLAA 5023 2r20	CGGAGTCTCCCACTTGCspGspC <sup>spT</sup>	0	1448,4	-	1751,6	1791,7	1766,6
71	HLAA 5024 2r20	GGCCGCCTCCCACTTGCspGspC <sup>spC</sup>	-14	1419,4	-	-	-	1737,6
72	HLAA 5271 1f20	AGTGGGAGACTCCGCCAspT <sup>spG</sup>	0	1255,3	1582,6	1558,5	-	1573,5
73	HLAA 5272 1f20	CAAGTGGGAGGCGGYYCCAspT <sup>spG</sup>	0	1255,3	1582,6	1558,5	-	1573,5
74	HLAA 5273 1f20	CAAGTGGGAGRCGGCCAspT <sup>spG</sup>	0	1255,3	1582,6	1558,5	-	1573,5
75	HLAA 5274 1f20	CAAGTGGGAGGCGGCCCTspT <sup>spG</sup>	0	1246,3	-	-	-	1564,5
76	HLAA 5275 1f20	CAAGTGGGAGGCGGCCGspT <sup>spT</sup>	0	1246,3	-	-	1589,6	-
77	HLAA 5276 1f20	CAAGTGGGAGGCGGCCGspT <sup>spC</sup>	0	1231,3	-	-	1574,5	-
78	HLAA 5277 1f20	CAAGTGGGAGGCGGCCMGspT <sup>spG</sup>	0	1271,3	1598,6	-	-	1589,5
79	HLAA 5278 1f20	CAAGTGGGAGGCRGCCGspT <sup>spG</sup>	0	1271,3	1598,6	-	-	1589,5
80	HLAA 5391 1f19	GCCCRTGAGGCGGAGCAspG <sup>spC</sup>	0	1138,1	1465,4	-	1481,4	1456,3
81	HLAA 5392 1f19	GYCCATGCGGCGGAGCAspG <sup>spC</sup>	0	1138,1	1465,4	-	1481,4	1456,3
82	HLAA 5393 1f19	GCCCGTCCGGCGGAGCAspG <sup>spC</sup>	0	1138,1	1465,4	-	1481,4	1456,3
83	HLAA 5394 1f19	GCCCATGTGGCGGAGCAspG <sup>spC</sup>	0	1138,1	1465,4	-	1481,4	1456,3
84	HLAA 5395 1f19	GTCCATGCGGCGGAGCAspG <sup>spT</sup>	0	1153,1	-	-	1496,4	1471,3
85	HLAA 5396 1f19	GCCCGTYGGGCGGAGCAspG <sup>spT</sup>	0	1153,1	-	-	1496,4	1471,3
86	HLAA 5397 1f19	GCCCATGAGGCGGAGCAspG <sup>spT</sup>	0	1153,1	-	-	1496,4	1471,3
87	HLAA 5398 1f19	GCCCWTGTGGCGGAGCAspG <sup>spT</sup>	0	1153,1	-	-	1496,4	1471,3
88	HLAA 5399 1f19	GCCMGTTGTGGCGGAGCAspG <sup>spT</sup>	0	1153,1	-	-	1496,4	1471,3
89	HLAA 5591 1r20	GCGGAGCCACTCCACGCAspC <sup>spT</sup>	0	1113,1	-	1416,3	-	-
90	HLAA 5592 1r20	GCGGAGCCCGTCCACGCAspC <sup>spT</sup>	0	1113,1	-	1416,3	-	-
91	HLAA 5593 1r20	GCGGAGCCACTCCACGCAspC <sup>spA</sup>	0	1122,1	-	-	1465,4	-
92	HLAA 5594 1r20	GCGGAGCCCGTCCACTCAspC <sup>spG</sup>	0	1138,1	-	-	-	1456,3
93	HLAA 5595 1r20	GCGGAGCCAGTCCACGCAspC <sup>spG</sup>	0	1138,1	-	-	-	1456,3
94	HLAA 5596 1r20	GCGGAGCCMGTCACGCAspC <sup>spG</sup>	0	1138,1	-	-	-	1456,3
95	HLAA 5597 1r20	GCGGAGCCACTCCACGCAspC <sup>spC</sup>	0	1098,1	1425,4	-	1441,4	-
96	HLAA 5598 1r20	GCGGAGCCCGTCCACGCAspC <sup>spC</sup>	0	1098,1	1425,4	-	1441,4	-
97	HLAA 5599 1r20	GCGGAGCCACTCCACGCAspG <sup>spG</sup>	0	1178,1	-	-	-	1496,3
98	HLAA 5711 2f20	TGGAGGGCCKGTGCGTGspGspA <sup>spG</sup>	0	1537,4	-	-	-	1855,6
99	HLAA 5712 2f20	TGGAGGGYAGTGCGTGspGspA <sup>spG</sup>	0	1537,4	-	-	-	1855,6
100	HLAA 5713 2f20	TGSAGGGCCGGTGCGTGspGspA <sup>spG</sup>	0	1537,4	-	-	-	1855,6
101	HLAA 5714 2f20	TGGATGSCACGTGCGTGspGspA <sup>spG</sup>	0	1537,4	-	-	-	1855,6
102	HLAA 5715 2f20	TGGAGGGCACSTGCGTGspGspA <sup>spG</sup>	0	1537,4	-	-	-	1855,6
103	HLAA 5716 2f20	TGGAGGGCACGTGMGTGspGspA <sup>spC</sup>	0	1497,4	-	-	1840,7	1815,6
104	HLAA 5717 2f20	TGGAGGGCYGGTGCGTGspGspA <sup>spC</sup>	0	1497,4	-	-	1840,7	1815,6

TABLE V

No	Name	Sequence	CT	Primer Masses	A	C	G	T
1	HLAB 971 2f20	CCCACTCCATGAGGCATspTspT <sup>Δ</sup> spC	0	1540,3	-	1843,7	1883,8	1858,7
2	HLAB 972 2f20	CCCACTYCATGAGGTATspTspT <sup>Δ</sup> spC	0	1540,3	-	1843,7	1883,8	1858,7
3	HLAB 2061 1f20	CGACGCCGCGAGTCMGAGspG <sup>Δ</sup> spA	-28	1150,1	1477,4	1453,3	-	1468,3
4	HLAB 2062 1f20	CGACGCCACGAGTCCGAGspG <sup>Δ</sup> spA	-28	1150,1	1477,4	1453,3	-	1468,3
5	HLAB 2063 1f20	CGACGCCGCGAGTCCRAGspA <sup>Δ</sup> spG	0	1178,1	1505,4	-	1521,4	-
6	HLAB 2064 1f20	CGACGCCRCGAGTCCGAGspA <sup>Δ</sup> spG	0	1178,1	1505,4	-	1521,4	-
7	HLAB 2221 1r19	GCCCCCTCTGCTCCACCspC <sup>Δ</sup> spA	0	1098,3	1425,4	-	1441,4	-
8	HLAB 2222 1r19	GCCCCCTCYTGCTCTATCspC <sup>Δ</sup> spA	0	1098,3	1425,4	-	1441,4	-
9	HLAB 2591 2f20	GGCCGGAGTATTGGGACspGspG <sup>Δ</sup> spG	0	1513,4	-	-	1856,7	-
10	HLAB 2592 2f20	GGCCGGAGTATTGGGACspGspA <sup>Δ</sup> spG	0	1497,4	-	-	1840,7	-
11	HLAB 2593 2f20	GGCCGGAGTATTGGGACspCspC <sup>Δ</sup> spG	-28	1405,4	-	-	1748,7	-
12	HLAB 2594 2f20	GGCCGGAGTATTGGGATspCspG <sup>Δ</sup> spG	0	1488,4	1815,7	-	1831,7	-
13	HLAB 2595 2f20	GGCCGGAGTTTTGGGACspCspG <sup>Δ</sup> spG	-28	1445,4	1772,7	-	1788,7	-
14	HLAB 2596 2f20	GGCCGGAGCATTGGGACspCspG <sup>Δ</sup> spG	-28	1445,4	1772,7	-	1788,7	-
15	HLAB 2597 2f20	GGCCGGGATATTGGGACspCspG <sup>Δ</sup> spG	-28	1445,4	1772,7	-	1788,7	-
16	HLAB 2598 2f20	GGCCRGAAATTGGGACspCspG <sup>Δ</sup> spG	-28	1445,4	1772,7	-	1788,7	-
17	HLAB 2599 2f20	GGCCGGMGTTATTGGGACspCspG <sup>Δ</sup> spG	-28	1445,4	1772,7	-	1788,7	-
18	HLAB 25910 2f20	GGCCTTAGTATTGGGACspCspG <sup>Δ</sup> spG	-28	1445,4	1772,7	-	1788,7	-
19	HLAB 2721 1f20	GGACSGGGAGACACGGAAspC <sup>Δ</sup> spA	0	1122,1	-	-	-	1440,3
20	HLAB 2722 1f20	GGACGRGGAGACACGGAAspC <sup>Δ</sup> spA	0	1122,1	-	-	-	1440,3
21	HLAB 2723 1f20	GGACCGGAACACACAGAAAspC <sup>Δ</sup> spT	0	1113,1	-	-	1456,4	-
22	HLAB 2724 1f20	GGACCGGAACACACAGACspC <sup>Δ</sup> spT	-14	1075,1	-	-	-	1393,3
23	HLAB 2725 1f20	GGACCGGGAGACACAGAAAspG <sup>Δ</sup> spT	0	1153,1	1480,4	-	-	-
24	HLAB 2726 1f20	GGACCGGGAGATACAGATspC <sup>Δ</sup> spT	0	1104,1	1431,4	1407,3	1447,4	1422,3
25	HLAB 2727 1f20	GGACCGGGASACACAGATspC <sup>Δ</sup> spT	0	1104,1	1431,4	1407,3	1447,4	1422,3
26	HLAB 2728 1f20	GGACCGGGACACACAGATspC <sup>Δ</sup> spT	0	1104,1	1431,4	1407,3	1447,4	1422,3
27	HLAB 2729 1f20	GGACCSGGAGACACAGATspC <sup>Δ</sup> spT	0	1104,1	1431,4	1407,3	1447,4	1422,3
28	HLAB 2921 2f19	CAAGACCAACACACAGspGspC <sup>Δ</sup> spT	0	1458,3	-	-	1801,6	-
29	HLAB 2922 2f19	CAAGSCCCAGGCACAGspGspC <sup>Δ</sup> spT	0	1458,3	-	-	1801,6	-
30	HLAB 2923 2f19	CAAGACCAACACACGspAspC <sup>Δ</sup> spT	-28	1414,3	-	-	1757,6	1732,5
31	HLAB 2924 2f19	GAAGGCCTCCGCGAGspAspC <sup>Δ</sup> spT	-28	1414,3	-	-	1757,6	1732,5
32	HLAB 2925 2f19	CAAGGCCMAGGCACAGspAspC <sup>Δ</sup> spT	-28	1414,3	-	-	1757,6	1732,5
33	HLAB 2926 2f19	CAAGSGCCAGGCACAGspAspC <sup>Δ</sup> spT	-28	1414,3	-	-	1757,6	1732,5
34	HLAB 2927 2f19	GAAGACCAACACACAGspAspC <sup>Δ</sup> spT	-28	1414,3	-	-	1757,6	1732,5
35	HLAB 3021 2f19	GCACAGACTGACCGAGspTspG <sup>Δ</sup> spG	0	1528,4	-	-	1871,7	-
36	HLAB 30211 2f19	ACACAGACTTACAGAGspAspG <sup>Δ</sup> spA	-28	1493,5	1820,8	-	1836,8	-
37	HLAB 3022 2f19	ACACAGACTTACCGAGspAspG <sup>Δ</sup> spG	0	1537,4	1864,7	-	-	-
38	HLAB 3023 2f19	RCACAGACTGACCGAGspAspG <sup>Δ</sup> spG	0	1537,4	1864,7	-	-	-
39	HLAB 3024 2f19	GCACAGACTGGCCGAGspTspG <sup>Δ</sup> spA	-28	1481,4	1811,7	-	1827,7	-
40	HLAB 3025 2f19	ACACAGACTTACCGAGspTspG <sup>Δ</sup> spA	-28	1481,4	1811,7	-	1827,7	-
41	HLAB 3026 2f19	RCACAGACTGACCGAGspTspG <sup>Δ</sup> spA	-28	1481,4	1811,7	-	1827,7	-
42	HLAB 3027 2f19	ACACAGGCTGACCGAGspAspG <sup>Δ</sup> spA	-28	1493,5	1820,8	-	1836,8	-
43	HLAB 3028 2f19	RCACAGACTGACCGAGspAspG <sup>Δ</sup> spA	-28	1493,5	1820,8	-	1836,8	-
44	HLAB 3029 2f19	GCRACAGACTTACCGAGspAspG <sup>Δ</sup> spA	-28	1493,5	1820,8	-	1836,8	-
45	HLAB 30210 2f19	ACACRGACTTACCGAGspAspG <sup>Δ</sup> spA	-28	1493,5	1820,8	-	1836,8	-
46	HLAB 3621 2f20	CGGGTCTCACACCTCCspAspC <sup>Δ</sup> spA	-28	1413,4	-	-	1756,7	-
47	HLAB 3622 2f20	CGGGTCTCACAYCATCCspAspG <sup>Δ</sup> spA	-14	1467,4	1794,7	1770,6	1810,7	1785,6
48	HLAB 3623 2f20	CGGKTCTCACACCTCCspAspG <sup>Δ</sup> spA	-14	1467,4	1794,7	1770,6	1810,7	1785,6
49	HLAB 3624 2f20	CGGGTCTCACACTTGGCspAspG <sup>Δ</sup> spA	-14	1467,4	1794,7	1770,6	1810,7	1785,6
50	HLAB 3625 2f20	CGGGTCTCACATCATCCspAspG <sup>Δ</sup> spG	-14	1483,4	-	-	-	1801,6
51	HLAB 3626 2f20	CGGGTCTCACACCTCCspAspG <sup>Δ</sup> spT	0	1472,4	-	-	1815,7	-
52	HLAB 3631 1r20	CCCASGTGCGAGCCGTACspA <sup>Δ</sup> spT	-28	1085,1	-	1388,3	1428,4	1403,3
53	HLAB 3632 1r20	CCCABGTGCGAGCCATACspA <sup>Δ</sup> spT	-28	1085,1	-	1388,3	1428,4	1403,3
54	HLAB 3633 1r20	CCCASGTGCGAGCCAAACspA <sup>Δ</sup> spT	-28	1085,1	-	1388,3	1428,4	1403,3

55	HLAB_3634_1r20	CCCACGTCGCAGCCAGACspA <sup>spT</sup>	-28	1085,1	-	1388,3	1428,4	1403,3
56	HLAB_3635_1r20	CCCACGTCGCAGCCGACspA <sup>spT</sup>	-28	1085,1	-	1388,3	1428,4	1403,3
57	HLAB_3636_1r20	CCCACGTCGCAGCCTTACspA <sup>spT</sup>	-28	1085,1	-	1388,3	1428,4	1403,3
58	HLAB_3637_1r20	CCCACGTCGCAGCCGTACspG <sup>spT</sup>	0	1129,1	-	1432,3	1472,4	1447,3
59	HLAB_3691_1f20	TCCGGCCCCAKGTCGCAGspC <sup>spC</sup>	0	1114,1	1441,4	-	1457,4	1432,3
60	HLAB_3692_1f20	TCCGGCCCCASGTCGCAGspC <sup>spC</sup>	0	1114,1	1441,4	-	1457,4	1432,3
55	HLAB_4121_2f20	GGCGCCTCCTCCGCGGGspTspA <sup>spC</sup>	-28	1444,4	-	1747,6	-	-
56	HLAB_4122_2f20	GGCGCCTCCTCCSCGGGspCspA <sup>spT</sup>	0	1472,4	1799,7	-	1815,7	-
57	HLAB_4123_2f20	GGCGCYTCTCCGCGGGspCspA <sup>spT</sup>	0	1472,4	1799,7	-	1815,7	-
58	HLAB_4124_2f20	GGCGTCTCCTCCGCGGTspTspA <sup>spT</sup>	0	1462,4	-	1765,6	-	-
59	HLAB_4125_2f20	GGCGCCTCCTCCGCGGGspTspA <sup>spT</sup>	-14	1473,4	-	1776,6	-	-
60	HLAB_4181_2f20	TCCTCCGCGGGTATGAAspCspA <sup>spG</sup>	0	1481,4	1808,7	-	-	-
61	HLAB_4182_2f20	TCCTCCACGGGTACCACspCspA <sup>spG</sup>	0	1457,4	-	-	-	1775,6
62	HLAB_4183_2f20	TCCTGCGCGGGTACCACspCspA <sup>spG</sup>	0	1457,4	-	-	-	1775,6
63	HLAB_4184_2f20	TCCTCCGCGGGTACCACspCspA <sup>spG</sup>	0	1457,4	-	-	-	1775,6
64	HLAB_4185_2f20	TCCTCTGCGGGTACCACspCspA <sup>spG</sup>	0	1457,4	-	-	-	1775,6
65	HLAB_4186_2f20	TCCTCCGCGGGTACCAGspCspA <sup>spG</sup>	0	1497,4	1824,7	1800,6	1840,7	1815,6
66	HLAB_4187_2f20	TMCTCCGCGGGTACCGGspCspA <sup>spG</sup>	0	1497,4	1824,7	1800,6	1840,7	1815,6
67	HLAB_4188_2f20	TCCTCCGCGGGTACCAGspCspG <sup>spG</sup>	0	1513,4	-	-	1856,7	-
68	HLAB_4191_2r20	AATCCTTGCCGTCGTAGspGspC <sup>spT</sup>	-14	1474,4	1801,7	-	-	-
69	HLAB_4192_2r20	AATCCTTGCCGTCGTAGspGspC <sup>spA</sup>	-28	1469,4	-	-	1812,7	-
70	HLAB_4193_2r20	AATCCTTGCCGTCGTAGspGspC <sup>spG</sup>	0	1513,4	1840,7	-	1856,7	1831,6
71	HLAB_4194_2r20	AATCTTTGCCGTCGTAGspGspC <sup>spG</sup>	0	1513,4	1840,7	-	1856,7	1831,6
72	HLAB_4195_2r20	AATCCTTGCCGTCGYAGspGspC <sup>spG</sup>	0	1513,4	1840,7	-	1856,7	1831,6
73	HLAB_4351n_1r20	TCMTTCAGGGCGATGTAAspT <sup>spC</sup>	-14	1201,3	-	1504,4	-	1519,4
74	HLAB_4352n_1r20	TCGTTTCAGGGCGATGTAAspT <sup>spT</sup>	0	1230,3	-	1533,5	-	-
75	HLAB_5271_1f20	CAAGTGGGAGGCGGCCCTspT <sup>spG</sup>	0	1246,3	-	-	-	1564,5
76	HLAB_5272_1f20	CAAGTKGGAGGCGGCCCGspT <sup>spG</sup>	0	1271,3	1598,6	1574,3	-	1589,5
77	HLAB_5391_1f20	GGCCCGTGYGGCGGAGCAspG <sup>spC</sup>	0	1138,1	-	-	1481,3	1456,3
78	HLAB_5392_1f20	GGCCCGTGTGCGGGAGCAspG <sup>spG</sup>	0	1178,1	1505,4	-	-	-
79	HLAB_5393_1f20	GGCCCGTGWGGCGGAGCAspG <sup>spG</sup>	0	1178,1	1505,4	-	-	-
80	HLAB_5394_1f20	GGCCCGTGAGGCGGAGCAspG <sup>spT</sup>	0	1153,1	-	-	1496,4	-
81	HLAB_5591_1r20	GCGGAGCGACTCCACGCAspC <sup>spT</sup>	0	1113,1	-	-	1456,4	-
82	HLAB_5592_1r20	GCGGAGCCACTCCACGCAspC <sup>spT</sup>	0	1113,1	-	-	1456,4	-
83	HLAB_5593_1r20	GCGGAGCCAATCCACGCAspC <sup>spT</sup>	0	1113,1	-	-	1456,4	-
84	HLAB_5594_1r20	GCGGAGCCACTCCACGCAspC <sup>spG</sup>	0	1152,1	-	-	-	1470,3
85	HLAB_5595_1r20	GCGGAGCGACTCCRCGCAspC <sup>spA</sup>	-14	1122,1	1449,1	1425,3	-	-
86	HLAB_5596_1r20	GCGGAGCSACTCCACGCAspC <sup>spA</sup>	-14	1122,1	1449,1	1425,3	-	-
87	HLAB_5597_1r20	GCGGAGCCCCTCCACGCAspC <sup>spA</sup>	-14	1122,1	1449,1	1425,3	-	-
88	HLAB_5711_1r20	CTCCAGGTAYCTGCGGAGspC <sup>spG</sup>	0	1154,1	1481,4	-	-	-
89	HLAB_5712_1r20	CTCCAGGTRTCTGCGGAGspC <sup>spC</sup>	0	1114,1	1441,4	1417,3	-	-
90	HLAB_583_1r19	ACCTGGAGAACGGGAAGspG <sup>spA</sup>	0	1178,1	1505,4	-	1521,4	-

TABLE VI

No	Name	Sequence	CT	Masses Primer	A	C	G	T
1	DRB1_1251_1r20	CATTGAAGAAATGACACTspC <sup>^</sup> spC	0	1098,1	-	1392,3	-	-
2	DRB1_1252_1r20	CGTTGAAGAAATGACACTspT <sup>^</sup> spA	0	1230,1	-	-	-	1548,5
3	DRB1_1253_1r20	CATTGAAGAAATGACATTspC <sup>^</sup> spA	0	1113,1	1440,4	1416,3	1456,4	1431,3
4	DRB1_1254_1r20	CATTGAAGAAWTAACACTspC <sup>^</sup> spA	0	1113,2	1440,4	1416,3	1456,4	1431,3
5	DRB1_1255_1r20	CRTTGAAGAAATGACACTspC <sup>^</sup> spA	0	1113,3	1440,4	1416,3	1456,4	1431,3
6	DRB1_1961_1f19	CATCTATAACCAAGAGGspA <sup>^</sup> spA	0	1162,1	-	-	-	1480,3
7	DRB1_1962_1f19	CTTCTATCACCAAGARGspA <sup>^</sup> spG	0	1178,1	1505,4	-	-	1496,3
8	DRB1_1963_1f19	CTTCTATAATCARGAGGspA <sup>^</sup> spG	0	1178,1	1505,4	-	-	1496,3
9	DRB1_1964_1f19	CGTCCATAACCAAGAGGspA <sup>^</sup> spG	0	1178,1	1505,4	-	-	1496,3
10	DRB1_1965_1f19	CATCTATAACCAAGAGGspA <sup>^</sup> spG	0	1178,1	1505,4	-	-	1496,3
11	DRB1_1966_1f19	CTTCCATAACCRGGAGGspA <sup>^</sup> spG	0	1178,1	1505,4	-	-	1496,3
12	DRB1_1967_1f19	CTTCGATAACCAGGAGGspA <sup>^</sup> spG	0	1178,1	1505,4	-	-	1496,3
13	DRB1_1968_1f19	CTTCTATAACCTGGAGGspA <sup>^</sup> spG	0	1178,1	1505,4	-	-	1496,3
14	DRB1_1971_1r20	CGTCGCTGTCTGAAGCGCAspG <sup>^</sup> spG	0	1178,1	1505,4	-	-	1496,3
15	DRB1_1972_1r20	CGTCGCTGTCTGAGCGCGspC <sup>^</sup> spG	0	1154,1	-	-	-	1472,3
16	DRB1_1973_1r20	CGTCGCTGTCTGAAGCGCAspA <sup>^</sup> spG	0	1162,1	-	-	-	1480,3
17	DRB1_1974_1r20	CGTCGCTGTCTGAAGYGCAspC <sup>^</sup> spG	-28	1110,1	1437,4	-	1453,4	1428,3
18	DRB1_1975_1r20	CGTCGCTGTCTGAASCGCAspC <sup>^</sup> spG	-28	1110,1	1437,4	-	1453,4	1428,3
19	DRB1_2271_1f20	CGACAGCGACGTGGGGGAspC <sup>^</sup> spT	0	1113,1	1440,4	-	-	-
20	DRB1_2272_1f20	CGACAGCGACGTGVGGGAspG <sup>^</sup> spT	0	1153,1	1480,4	-	-	1471,3
21	DRB1_2611_1r20	TTCTGGCTGTTCCAGTACspT <sup>^</sup> spG	0	1231,2	-	-	1574,5	-
22	DRB1_2612_1r20	TTCTGGCTGTTCCAGTACspC <sup>^</sup> spC	0	1074,1	-	1377,3	-	-
23	DRB1_2613_1r20	TTCTGGCTGTTCCAGTAGspT <sup>^</sup> spC	0	1231,2	-	1534,4	-	-
24	DRB1_2614_1r20	TTCTGGCTGTTCCAGTRCspT <sup>^</sup> spC	-14	1177,2	1504,5	1480,4	1520,5	-
25	DRB1_2615_1r20	TTCTGGCTGTTCCAGGACspT <sup>^</sup> spC	-14	1177,2	1504,5	1480,4	1520,5	-
26	DRB1_2861_1f19	CTGGAACAGCCAGAAGAspA <sup>^</sup> spC	-28	1122,1	1449,4	-	-	-
27	DRB1_2862_1f19	CTGGAACAGCCRGAAGGspA <sup>^</sup> spC	0	1138,1	1465,4	1441,3	-	1456,3
28	DRB1_2991_1f20	GAAGGACHTCCTGGAGCAspG <sup>^</sup> spG	0	1178,1	-	1481,3	-	-
29	DRB1_2992_1f20	GAAGGACATCCTGGGAGAspC <sup>^</sup> spA	-14	1108,1	1435,1	-	1451,4	-
30	DRB1_2993_1f20	GAAGGACATCCTGGARGAspC <sup>^</sup> spA	-14	1108,1	1435,1	-	1452,4	-
31	DRB1_2994_1f20	GAAGGACYTCCTGGAAGAspC <sup>^</sup> spA	-14	1108,1	1435,1	-	1453,4	-
32	DRB1_2995_1f20	GAAGGACATCCTGGAGCAspG <sup>^</sup> spA	0	1162,1	1489,4	-	1505,4	-
33	DRB1_2996_1f20	GAAGGACHTCCTGGAGCGspG <sup>^</sup> spA	0	1178,1	-	-	1521,4	-
34	DRB1_2997_1f20	GAAGGACHTCCTGGAAGAspC <sup>^</sup> spG	0	1138,1	1465,4	-	-	-
35	DRB1_3081_1r20	GTCTGCAATAGGTGTCCAspC <sup>^</sup> spG	0	1138,1	-	1441,3	-	-
36	DRB1_3082_1r20	GTCTGCARTAGGCGTCCAspC <sup>^</sup> spC	-14	1084,1	1411,4	1387,3	1427,4	1402,3
37	DRB1_3083_1r20	GTCTGCAGTAATTGTCCAspC <sup>^</sup> spC	-14	1084,1	1411,4	1387,3	1427,4	1402,3
38	DRB1_3084_1r20	GTCTGCACACGGTGTCCAspC <sup>^</sup> spC	-14	1084,1	1411,4	1387,3	1427,4	1402,3
39	DRB1_3085_1r20	GTCTGCAGTAGGTGTCCAspC <sup>^</sup> spC	-14	1084,1	1411,4	1387,3	1427,4	1402,3
40	DRB1_3086_1r20	GTCTGCAATAGGTGTCCAspC <sup>^</sup> spC	-14	1084,1	1411,4	1387,3	1427,4	1402,3
41	DRB1_341_1f19	TGCAGACACAACACTACSGspG <sup>^</sup> spG	0	1194,1	-	1497,3	-	1512,3
42	DRB1_3451_1r20	CGCTGCACTGTGAATCTCspT <sup>^</sup> spC	0	1191,3	1518,5	1494,4	-	-
43	DRB1_3452_1r20	CTCTGCACTGTGAAGCTCspT <sup>^</sup> spC	0	1191,3	1518,5	1494,4	-	-
44	DRB1_3453_1r20	CGCTGCACYGTGAAGCTCspT <sup>^</sup> spC	0	1191,3	1518,5	1494,4	-	-

The resolution achievable by 19 markers each for HLA-A and HLA-B and the ten markers for HLA-DRB1 are listed in Tables VII to IX below.

**TABLE VII**

Frequent Alleles of HLA-A	Group of frequent Alleles with same four-digit type	Rare Alleles with same Mini-Haplotype Profile	Resolution (in %)
A*0101	A*010101, A*010102	A*0103, A*0104N, A*0109	98,3
A*0201	A*02010101, A*02010102L, A*020103, A*020104, A*020108, A*020109	A*0204, A*0209, A*0225, A*0231, A*0232N, A*0242, A*0243N, A*0253N, A*0258, A*0260, A*0264, A*0266, A*0267	93,4
	A*020102		100
	A*020105		100
	A*020106		100
	A*020107		100
A*0301	A*03010101, A*03010102N	A*0303N, A*0304, A*0305, A*0306, A*0311N	97,6
	A*030102		100
	A*030103		100
A*2301	A*2301	A*2306, A*2307N, A*2308N	98,6
A*2402	A*24020101, A*24020102L, A*240202, A*240203, A*240204	A*2404, A*2409N, A*2411N, A*2426, A*2427, A*2432, A*2435, A*2436N, A*2437, A*2439	94,5
A*2902	A*290201	A*29010101, A*29010102N, A*2906, A*2908N	98,3
	A*290202		100
A*3001	A*3001		100
A*3002	A*3002		100

5

Capture: Alleles in a same field have the same mini-haplotype profile; grey highlighted are all alleles with identical sequences over exons 2 and 3.

10

15

TABLE VIII

Frequent Alleles of HLA-B	Groups of frequent Alleles with same four-digit type	Rare Alleles with same Mini-Haplotype Profile	Resolution (in %)
B*0702	B*070201, B*070202, B*070203, B*070204	B*0703, B*0721, B*0722, B*0723, B*0730, B*0733, B*0735	98,0
B*0801	B*0801	B*0808N, B*0818, B*0819N	99,3
B*1302	B*1302	B*1308	99,6
B*1501	B*15010101, B*15010102N, B*150103, B*150104	B*1528, B*1533, B*1534, B*1560, B*1575, B*1578, B*1579N, B*1581, B*1582	97,6
	B*150102		100
B*1801	B*180101, B*180102	B*1805, B*1817N	99,3
B*3501	B*350101, B*350102	B*3507, B*3540N, B*3541, B*3542, B*5305	98,7
B*3503	B*3503	B*3536	99,6
B*4001	B*400101, B*400102	B*4011, B*401401, B*401402, B*401403, B*4022N	98,7
	B*400103		100
	B*400104	B*4004	99,6
B*4402	B*44020101, B*44020102S, B*440202, B*440203	B*4411, B*4419N, B*4422, B*4423N, B*4427, B*4433, B*4434, B*4435	97,8
B*4403	B*440301	B*4413, B*4426, B*4429, B*4430, B*4432, B*4436, B*4437, B*4438, B*4439	98,2
	B*440302	B*4407	99,6
B*5101	B*510101, B*510102, B*510105	B*5111N, B*5112, B*5114, B*5118, B*5126, B*5127N, B*5128, B*5130, B*5132, B*5133	97,6
	B*510103		100
	B*510104	B*5124	99,6
B*5701	B*570101	B*5706, B*5708	99,5
	B*570102		100

Capture: Alleles in a same field have the same mini-haplotype profile; grey highlighted are all alleles with identical sequences over exons 2 and 3.

**TABLE IX**

Frequent Alleles of HLA-DRB1*	Groups of frequent Alleles with same four-digit type	Rare Alleles with same Mini-Haplotype Profile	Resolution (in %)
DRB1*0101	DRB1*010101	DRB1*0105, DRB1*0107, DRB1*0111	98,9
	DRB1*010102		100
DRB1*0301	DRB1*030101, DRB1*030102	DRB1*0307, DRB1*0312, DRB1*0313, DRB1*0315, DRB1*0316, DRB1*0318, DRB1*0322, DRB1*0323	97,2
DRB1*0401	DRB1*040101, DRB1*040102	DRB1*0409, DRB1*0426, DRB1*0433	98,6
DRB1*0701	DRB1*070101, DRB1*070102	DRB1*0703, DRB1*0704, DRB1*0705, DRB1*0707	98,3
DRB1*1101	DRB1*110101, DRB1*110102, DRB1*110103, DRB1*110104, DRB1*110105	DRB1*112701, DRB1*112702, DRB1*1130, DRB1*1139	97,5
DRB1*1104	DRB1*110401, DRB1*110402	DRB1*1134, DRB1*1146	98,9
DRB1*1302	DRB1*130201, DRB1*130202	DRB1*1331, DRB1*1339, DRB1*1341	98,6
DRB1*1501	DRB1*150101, DRB1*150103, DRB1*150105	DRB1*1503, DRB1*1506, DRB1*1509, DRB1*1513	98,0
	DRB1*150102		100
	DRB1*150104	DRB1*1512	99,4

Capture: Alleles in a same field have the same mini-haplotype profile; grey highlighted are all alleles with identical sequences over exon 2 (base 101 to 356)

- 5 The complete list of HLA alleles and sub-groups generated by the most informative mini-haplotyping markers (ten each for HLA-A, HLA-B and HLA-DRB1) are listed in Tables X to XII below.

### TABLE X

[illegible]



	A*1109	T	C	T	A	C	C	G	G	A	G	C	A	C	C	A	G	G	G	A	G	T	T	G	T	A	C	A	G	G	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*1112	T	C	T	A	C	C	G	G	A	G	C	A	C	C	A	G	G	G	A	G	T	T	G	T	A	C	A	G	G	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*1115	T	C	T	A	C	C	G	G	A	G	C	A	C	C	A	G	G	G	A	G	T	T	G	T	A	C	A	G	G	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*1113	T	C	T	A	C	C	G	G	A	G	C	A	C	C	A	G	G	G	A	G	T	T	G	T	A	C	A	G	G	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*1105	T	C	T	A	C	C	G	G	A	G	C	A	C	C	A	G	G	G	A	G	T	T	G	T	A	C	A	G	G	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*1114	T	C	T	A	C	C	G	G	A	G	C	A	C	C	A	G	G	G	A	G	T	T	G	T	A	C	A	G	G	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*1110	T	C	T	A	C	C	G	G	A	G	C	A	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6809	T	C	T	A	C	C	G	G	A	G	C	A	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
5	A*1111	T	C	T	A	C	C	G	G	A	G	C	A	C	C	A	G	G	G	A	G	T	T	G	T	A	C	T	G	C	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*1108	T	C	T	A	C	C	G	G	A	G	C	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	A	G	G	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6805	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	C	G	G	G	C	C	T	G	T	G		
	A*6820	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	C	G	G	G	C	C	T	G	T	G		
	A*680301	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*680302	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6804	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	T	T	G	G	G	G	C	C	T	G	T	G		
	A*680101	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*680102	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*680103	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6807	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6811N	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6812	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6817	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6819	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6821	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6822	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6823	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6824	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6816	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6813	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
15	A*6810	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	G	A	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6814	T	C	T	A	C	C	G	G	A	G	T	G	C	C	A	G	G	G	A	G	T	T	G	T	A	G	A	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*6808	T	C	T	A	C	C	G	G	A	G	C	T	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*0312	T	C	T	A	C	C	G	G	A	G	T	T	C	C	A	G	G	G	A	G	T	T	G	T	A	C	A	G	G	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*3402	T	C	T	A	C	C	G	G	A	G	T	T	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*3403	T	C	T	A	C	C	G	G	A	G	T	T	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*3404	T	C	T	A	C	C	G	G	A	G	T	T	C	C	A	G	G	G	A	G	T	T	G	T	A	C	G	G	A	A	C	T	G	A	G	G	C	C	T	G	T	G		
20	A*3102	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	A	G	G	A	A	T	T	G	A	G	G	C	C	A	G	T	G	
	A*3107	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	A	G	G	A	A	C	T	G	A	G	G	C	C	A	G	T	G	
	A*3108	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	G	A	G	G	A	A	C	T	G	G	G	G	C	C	A	G	T	G	
	A*3105	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	C	G	T	T	G	T	A	C	A	G	G	A	A	T	T	G	A	G	G	C	C	T	G	T	G
	A*310102	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	A	G	G	A	A	T	T	G	A	G	G	C	C	T	G	T	G	
	A*3109	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	A	G	G	A	A	T	T	G	A	G	G	C	C	T	G	T	G	
	A*3301	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	T	T	G	G	G	G	C	C	T	G	T	G		
	A*330301	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	T	T	G	G	G	G	C	C	T	G	T	G		
	A*330302	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	T	T	G	G	G	G	C	C	T	G	T	G		
	A*3304	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	T	T	G	G	G	G	C	C	T	G	T	G		
	A*3305	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	T	T	G	G	G	G	C	C	T	G	T	G		
	A*3306	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	T	T	G	G	G	G	C	C	T	G	T	G		
	A*3307	T	C	A	C	C	C	A	G	A	G	T	T	C	C	A	C	G	G	A	G	T	T	G	T	A	C	G	G	A	A	T	T	G	G	G	G	C	C	T	G	T	G		
	A*2905	T	C	A	C	C	C	G	G	A	G	C	A	C	C	A	G	G	G	A	G	T	T	G	T	A	C	T	G	C	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*2903	T	C	A	C	C	C	G	G	A	G	T	T	C	C	A	G	G	G	A	C	G	T	T	G	A	C	T	G	C	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*290202	T	C	A	C	C	C	G	G	A	G	T	T	C	C	A	G	G	G	A	G	T	T	G	T	A	C	T	G	C	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*290101	T	C	A	C	C	C	G	G	A	G	T	T	C	C	A	G	G	G	A	G	T	T	G	T	A	C	T	G	C	A	C	T	G	G	G	G	C	C	T	G	T	G		
	A*290101	T	C	A	C	C	C	G	G	A	G	T	T	C	C	A	G	G	G	A	G	T	T	G	T	A	C	T	G	C	A	C	T	G	G	G	G	C	C	T	G	T	G		
</																																													

A*3104	T	C	A	C	C	G	G	A	G	T	T	C	C	A	G	G	A	G	T	T	T	G	T	A	C	A	G	G	A	T	T	G	A	G	G	C	C	T	G	T	G				
A*3106	T	C	A	C	C	G	G	A	G	T	T	C	C	A	G	G	A	G	T	T	T	G	T	A	C	A	G	G	A	T	T	G	A	G	G	C	C	T	G	T	G				
A*3004	T	T	C	T	T	G	A	A	A	G	T	G	C	C	A	C	G	G	A	G	T	T	T	G	T	A	C	A	G	G	A	C	T	G	A	G	G	C	C	T	G	T	G		
A*3006	T	T	C	T	C	T	G	A	A	A	G	T	G	C	C	A	C	G	G	A	G	T	T	T	G	T	A	C	A	G	G	A	C	T	G	A	G	G	C	C	T	G	T	G	
A*3010	T	T	C	T	C	T	G	A	A	A	G	T	T	C	C	A	C	G	G	A	G	T	T	T	G	C	A	C	A	G	G	A	C	T	G	A	G	G	C	C	T	G	T	G	
A*3007	T	T	C	T	C	T	G	A	A	A	G	T	T	C	C	A	C	G	G	A	G	T	T	T	G	T	A	G	A	G	G	A	C	T	G	A	G	G	C	C	A	G	T	G	
A*3002	T	T	C	T	C	T	G	A	A	A	G	T	T	C	C	A	C	G	G	A	G	T	T	T	T	G	T	A	C	A	G	G	A	C	T	G	A	G	G	C	C	T	G	T	G
A*3009	T	T	C	T	C	T	G	A	A	A	G	T	T	C	C	A	C	G	G	A	G	T	T	T	T	G	T	A	C	A	G	G	A	C	T	G	A	G	G	C	C	T	G	T	G
A*3012	T	T	C	T	C	T	G	A	A	A	G	T	T	C	C	A	C	G	G	A	G	T	T	T	T	G	T	A	C	A	G	G	A	C	T	G	A	G	G	C	C	T	G	T	G
A*3003	T	T	C	T	C	T	G	A	A	A	G	T	T	C	C	A	C	G	G	A	G	T	T	T	T	G	T	A	C	A	G	G	A	C	T	G	G	G	C	C	T	G	T	G	
A*3001	T	T	C	T	C	T	G	A	A	A	G	T	T	C	C	A	C	G	G	A	G	T	T	T	T	G	T	A	C	A	G	G	A	C	T	G	A	G	G	C	C	T	G	T	G
A*3011	T	T	C	T	C	T	G	A	A	A	G	T	T	C	C	A	C	G	G	A	G	T	T	T	T	G	T	A	C	A	G	G	A	C	T	G	G	G	C	C	T	G	T	G	
A*2414	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	G	G	A	C	T	G	G	G	C	C	A	G	T	G		
A*2415	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	G	G	A	C	T	G	G	G	C	C	A	G	T	G		
A*2428	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	G	G	A	C	T	G	G	G	C	C	A	G	T	G		
A*2430	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	G	G	A	C	T	G	G	G	C	C	A	G	T	G		
A*2408	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	G	G	A	G	G	A	C	T	G	G	G	C	C	A	G	T	G		
A*2431	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2420	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2429	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*240201	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*240201	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*240202	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*240203	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*240204	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2404	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2409N	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2411N	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2426	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2427	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2432	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2435	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2436N	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2437	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2438	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2439	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2421	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2425	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2405	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2434	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	T	T	G	G	G	C	C	A	G	T	G					
A*2407	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2419	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2423	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	C	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G						
A*240301	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G						
A*240302	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G						
A*2433	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G						
A*2410	T	T	C	T	C	C	C	A	C	A	G	C	A	C	C	A	C	G	G	A	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G						
A*2302	T	T	C	T	C	C	C	A	C	A	G	T	G	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2406	T	T	C	T	C	C	C	A	C	A	G	T	G	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2422	T	T	C	T	C	C	C	A	C	A	G	T	G	C	C	A	C	G	G	A	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G						
A*2301	T	T	C	T	C	C	C	A	C	A	G	T	T	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2306	T	T	C	T	C	C	C	A	C	A	G	T	T	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2307N	T	T	C	T	C	C	C	A	C	A	G	T	T	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2308N	T	T	C	T	C	C	C	A	C	A	G	T	T	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2305	T	T	C	T	C	C	C	A	C	A	G	T	T	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2310	T	T	C	T	C	C	C	A	C	A	G	T	T	C	C	A	C	G	G	A	C	G	T	T	T	G	A	G	A	C	T	G	G	G	C	C	A	G	T	G					
A*2413	T	T	C	T	C	C	C	A	C	A	G	T	T	C	C	A	C	G	G	A																									

5	A*2616	TCTC	CCAG	AGTG	CCAC	GGAG	TGTA	CGGA	ACTG	GGGC	CTGT	G
	A*3207	TCTC	CCAG	AGTT	CCAC	GGAG	TGTA	CAGG	ACTG	GGGC	CTGT	G
	A*0102	TCTC	CCGG	AGCG	CCAC	GGAC	TGTA	CAGG	ACTG	GGGC	CTAT	G
	A*2417	TCTC	CCGG	AGCA	CCAC	GGAC	TGTT	GAGG	ACTG	GGGC	CAGT	G
	A*0252	TCTT	TGAA	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
10	A*0219	TCTT	CCAC	AGCA	CCAC	GGAC	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0237	TCTT	CCAC	AGCA	CCAC	GGAC	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0213	TCTT	CCAC	AGCA	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0238	TCTT	CCAC	AGCA	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0227	TCTT	CCAC	AGCA	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
15	A*0212	TCTT	CCAC	AGCA	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0249	TCTT	CCAC	AGCG	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0250	TCTT	CCAC	AGTG	CCAC	GGAG	TGTA	GGGG	ACTC	GGGC	CAGT	G
	A*0203	TCTT	CCAC	AGTG	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0202	TCTT	CCAC	AGTG	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
20	A*0222	TCTT	CCAC	AGTG	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0263	TCTT	CCAC	AGTG	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0247	TCTT	CCAC	AGTG	CCAC	GGAG	TGTA	GGGG	AGTC	GGTC	CAGT	G
	A*0236	TCTT	CCAC	AGTT	CCAC	GGAC	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0245	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	CAGG	ACTC	GGTC	CAGT	G
25	A*0246	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GAGG	ACTC	GGTC	CAGT	G
	A*0248	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GAGG	ACTG	GGGC	CAGT	G
	A*0255	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	CGGA	ACTC	GGGC	CTGT	G
	A*020105	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0226	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
30	A*0259	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*020106	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*020101	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*020101	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*020102	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
35	A*020103	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*020104	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*020108	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*020109	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0204	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
40	A*0209	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0225	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0231	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0232N	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0242	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
45	A*0243N	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0253N	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0258	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0260	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0264	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
50	A*0266	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0267	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*020107	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0216	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0229	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
55	A*0230	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0240	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*0224	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*022002	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
	A*022001	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ACTC	GGTC	CAGT	G
60	A*0211	TCTT	CCAC	AGTT	CCAC	GGAG	TGTA	GGGG	ATTG	GGTC	CAGT	G
	A*0233	TCTT	CCAC	AGTT	CCAC	GGAG	TGTC	GGGG	ACTC	GGTC	CAGT	G
	A*0207	TCTT	CCAC	AGTT	CCAC	GGAG	TGTC	GGGG	ACTC	GGTC	CAGT	G
	A*0215N	TCTT	CCAC	AGTT	CCAC	GGAG	TGTC	GGGG	ACTC	GGTC	CAGT	G
	A*0218	TCTT	CCAC	AGTT	CCAC	GGAG	TGTC	GGGG	ACTC	GGTC	CAGT	G
65	A*021701	TCTT	CCAC	AGTT	CCAC	GGAG	TGTT	GGGG	ACTC	GGTC	CAGT	G

[illegible]



B-1551	A	G	G	A	C	C	G	C	C	C	T	G	T	G	C	A	T	G	T	C	T	G	C	A	G	A	G	G	A	G	A	T	G	G	A	A	G	T	A	C			
B-4417	A	G	G	A	C	C	G	C	C	C	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
B-4418	A	G	G	A	T	A	G	C	C	C	T	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
B-4415	A	G	G	A	T	A	G	C	C	C	C	T	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C	
B-4501	A	G	G	A	T	A	G	C	C	C	C	T	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C	
B-4503	A	G	G	A	T	A	G	C	C	C	C	T	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C	
B-4504	A	G	G	A	T	A	G	C	C	C	C	T	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C	
B-4505	A	G	G	A	T	A	G	C	C	C	C	T	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C	
B-4506	A	G	G	A	T	A	G	C	C	C	C	T	G	T	G	T	A	T	A	T	C	T	A	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C	
B-8201	A	G	G	A	T	A	G	C	C	C	C	T	G	T	G	C	A	T	A	T	C	T	A	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	T	C	
B-8202	A	G	G	A	T	A	G	C	C	C	C	T	G	T	G	C	A	T	A	T	C	T	A	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	T	C	
B-3538	A	G	G	A	T	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	C
B-3701	A	G	G	A	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	C	T
B-3703N	A	G	G	A	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	C	T
B-3704	A	G	G	A	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	C	T
B-3705	A	G	G	A	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	C	T
B-4502	A	G	G	A	T	C	G	C	C	C	C	T	G	T	G	T	A	T	A	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	A	G	T	A	C
B-0731	A	G	C	G	A	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	A	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
B-0806	A	G	C	G	A																																						

B-0710	A	G	C	G	A	C	G	C	C	G	A	G	T	G	C	A	T	G	T	C	T	G	C	A	G	A	G	G	A	G	A	G	A	G	T	A	C							
B-0708	A	G	C	G	A	C	G	C	C	G	A	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	A	G	A	G	T	A	C							
B-0714	A	G	C	G	A	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	G	G	A	G	A	G	A	G	T	A	C							
B-0728	A	G	C	G	A	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	G	G	A	G	A	G	A	G	T	A	C							
B-0709	A	G	C	G	C	C	G	C	C	G	A	G	T	G	C	A	T	G	T	C	T	A	C	A	G	A	G	G	A	G	A	G	A	G	T	A	C							
B-0711	A	G	C	G	C	C	G	C	C	G	A	G	T	G	C	A	T	G	T	C	T	A	C	A	G	A	G	G	A	A	C	A	T	G	G	A	G	T	A	C				
B-1547	A	G	C	G	C	C	G	C	C	G	A	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	G	A	T	G	G	A	G	T	A	C			
B-1549	A	G	C	G	C	C	G	C	C	G	A	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	G	A	T	G	G	A	G	T	A	C			
B-0717	A	G	C	G	C	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	G	G	A	G	A	C	A	T	G	G	A	G	T	A	C			
B-1565	A	G	C	G	C	C	G	C	C	C	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	G	A	T	G	G	G	A	T	A	C			
B-1813	A	G	C	G	C	C	G	C	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	G	A	T	G	G	G	A	C	T	A	C	
B-3508	A	G	C	G	C	C	G	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	C	A	G	A	G	G	A	G	A	G	A	T	G	G	G	A	C	G	T	A	C
B-3545	A	G	C	G	C	C	G	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	C	A	G	A	G	G	A	G	A	G	A	T	G	G	G	A	C	G	T	A	C
B-5509	A	G	C	G	T	A	G	C	C	C	G	A	G	T	G	C	A	T	A	T	C	T	A	C	A	G	A	C	G	A	G	A	G	A	T	G	G	A	G	A	T	A	C	
B-1403	A	G	C	G	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	T	C	T	G	C	A	G	T	G	G	A	G	A	G	A	T	G	G	A	G	A	T	A	C	
B-3908	A	G	C	G	T	C	G	C	C	C	A	C	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	G	G	A	G	A	G	A	T	G	G	A	G	A	T	A	C	
B-3911	A	G	C	G	T	C	G	C	C	C	A	C	G	T	G	C	A	T	A	T	C	T	G	C	A	G	A	G	G	A	G	A	G	A	T	G	G	A	G	A	T	A	C	
B-1530	A	G	T	G	A	C	G	C	C	C	T	G	T	G	C	A	T	A	T	C	T	C	C	C	A	G	A	C	G	A	G	A	G	A	T	G	G	G	A	T	G	T	A	C
B-1545	A	G	T	G	A	C	G	C	C	C	T	G	T	G	C	A	T	G	T	C	T	C	C	C	A	G</																		

5

10

15

20

25

30

B-1578	A	G	T	G	C	C	G	C	C	C	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	A	C	C			
B-1579N	A	G	T	G	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	A	C	C		
B-1581	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	A	C	C	
B-1582	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	A	C	C	
B-1527	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	T	T		
B-1532	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	C	T		
B-1557	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	T	G	A	A	G	A	T	G	G	G	A	T	G	T	A	C
B-1566	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	G	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	A	C		
B-1508	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	A	C		
B-1556	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	A	C		
B-351401	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	C	G	T	A	T		
B-351402	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	C	G	T	A	T		
B-3543	A	G	T	G	C	C	G	C	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C		
B-1573	A	G	T	G	T	A	G	C	C	C	C	C	T	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	A	C		
B-1558	A	G	T	G	T	C	G	C	C	C	C	C	T	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	G	A	T	G	T	A	C		
B-3918	A	G	T	G	T	C	G	C	C	C	C	C	A	C	G	T	G	C	A	T	A	T	C	T	G	C	A	G	A	G	G	A	G	G	A	T	G	G	A	G	A	G	T	A	C	
B-0734	A	G	C	T	A	C	G	C	C	C	C	C	A	C	G	T	G	C	A	T	A	T	C	T	A	C	A	G	A	G	G	A	G	C	A	T	G	G	A	G	A	G	T	A	C	
B-5504																																														



5	B-5114	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-5118	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-5119	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-5123	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-5124	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-5126	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-5127N	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-5128	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-5130	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-5132	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
B-5133	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
10	B-7801	A	G	C	T	A	C	G	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
	B-780201	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-780202	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-7804	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-3502	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-3504	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-350901	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-350902	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-3512	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
	B-3522	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
B-5104	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
B-5106	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
B-5117	A	G	C	T	A	C	G	C	C	C	T	G	T	T	T	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T		
B-15170101	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	C	C	A	C	A	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T		
B-15170102	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	C	C	A	C	A	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T		
B-1510	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
B-1537	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
B-3534	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
B-3539	A	G	C	T	A	C	G	C	C	C	T	G	T	G	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
15	B-440301	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
	B-4413	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
	B-4426	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
	B-4429	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
	B-4430	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
	B-4432	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
	B-4436	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
	B-4437	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
	B-4438	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
	B-4439	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
B-440302	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
B-4407	A	G	C	T	A	C	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T	
B-8101	A	G	C	T	A	C	G	C	C	C	G	A	G	T	G	G	C	A	T	A	T	C	T	A	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
B-8102	A	G	C	T	A	C	G	C	C	C	G	A	G	T	G	G	C	A	T	A	T	C	T	A	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
B-40060101	A	G	C	T	A	C	G	C	C	C	G	A	G	T	G	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
B-40060102	A	G	C	T	A	C	G	C	C	C	G	A	G	T	G	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
B-4019	A	G	C	T	A	C	G	C	C	C	G	A	G	T	G	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
B-4037	A	G	C	T	A	C	G	C	C	C	G	A	G	T	G	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
B-4047	A	G	C	T	A	C	G	C	C	C	G	A	G	T	G	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
B-4431	A	G	C	T	A	C	G	C	C	C	G	A	G	T	G	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	A	G	A	T	G	G	G	A	C	G	G	T	A	T
B-4002	A	G	C	T	A	C	G	C	C	C	G	A	G	T	G	G	C	A	T	A	T	C	T	C	C	A	G	A																			

	B-4004	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	G	G	A	G	A	G	G	A	T	G	G	A	A	G	T	A	T	T
	B-4012	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	A	A	G	T	A	T	T		
	B-4046	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	A	A	G	T	A	T	T		
	B-4803	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	G	G	A	G	G	A	T	G	G	A	A	G	T	A	T	T		
	B-4030	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	G	T	G	A	G	G	A	T	G	G	A	A	G	T	A	T	T		
	B-4034	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	G	T	G	A	G	G	A	T	G	G	A	A	G	T	A	T	T		
	B-2720	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	G	C	A	G	A	A	G	A	G	T	A	T	G	G	A	A	G	T	A	T	T		
	B-2707	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	G	C	A	G	A	G	G	A	G	C	A	T	G	G	A	A	G	T	A	T	T		
5	B-2711	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	G	C	A	G	A	G	G	A	G	C	A	T	G	G	A	A	G	T	A	T	T		
	B-2724	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	G	C	A	G	A	G	G	A	G	C	A	T	G	G	A	A	G	T	A	T	T		
	B-5110	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	T	G	G	A	A	G	T	A	T	T			
	B-5116	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	T	G	G	A	A	G	T	A	T	T			
	B-5131	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	T	G	G	A	A	G	T	A	T	T			
	B-5134	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	T	C	A	G	A	C	G	A	G	A	T	G	G	A	A	G	T	A	T	T			
	B-3531	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	T	C	A	G	A	G	G	A	G	C	A	T	G	G	A	A	G	T	A	T	T		
	B-4007	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	A	G	T	A	T	T		
	B-4008	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	T	C	A	G	A	G	G	A	G	C	A	T	G	G	A	A	G	T	A	T	T		
10	B-4013	A	G	C	T	A	C	G	C	C	G	A	G	T	G	C	A	T	A	T	C	T	T	C	A	G	A	G	G	A	A	C	A	T	G	G	A	A							

	B-180101	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
	B-180102	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
	B-1803	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
	B-1804	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
	B-1805	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
	B-1811	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
	B-1812	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
	B-1815	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
	B-1817N	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
5	B-1808	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	G	C			
	B-1818	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	C	T			
	B-1806	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	T	G	A	G	G	A	T	G	G	G	A	C		
	B-3907	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	G	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
	B-1544	A	G	C	T	C	C	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	G	C	A	G	A	G	G	A	T	G	G	G	A	T	G	T	A	T			
	B-1807	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	C			
	B-3535	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	T			
10	B-5805	A	G	C	T	C	C	G	C	C	A	C	G	T	G	C	A	T	G	A	C	A	T	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	T
	B-5609	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	A	C	A	G	A	G	G	A	T	G	G	G	A	C	G	T	A	T			
	B-1304	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	T	G	T	A	T
	B-5309	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	T
	B-1503	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-1561	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-1564	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-1574	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-1505	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
15	B-1539	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-1562	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	T	
	B-4802	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	T	
	B-1520	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	T	
	B-3520	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	T	
	B-3528	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	T	
	B-1531	A	G	C	T	C	C	G	C	C	G	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-1555	A	G	C	T	C	C	G	C	C	G	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-1513	A	G	C	T	C	C	G	C	C	G	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	T
20	B-1536	A	G	C	T	C	C	G	C	C	G	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	T
	B-1502	A	G	C	T	C	C	G	C	C	G	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	T	
	B-1525	A	G	C	T	C	C	G	C	C	G	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	T	
	B-1506	A	G	C	T	C	C	G	C	C	G	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	T	T	
	B-1523	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	G	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	C
	B-1518	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	G	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-1572	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	G	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-3526	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	G	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	T	
	B-1521	A	G	C	T	C	C	G	C	C	G	T	G	T	G	T	A	T	G	T	C	T	G	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	T	
25	B-1516	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	T
	B-1567	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	T
	B-3537	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	T	
	B-1529	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-3505	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-3516	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-3517	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-3530	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-3532	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	A	T	G	G	G	A	C	G	T	A	C	
	B-3527	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	T
30	B-5301	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	T
	B-5302	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C	G	T	A	T
	B-5306	A	G	C	T	C	C	G	C	C	G	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	A	G	A	T	G	G	G	A	C				

	B-3525	A	G	C	T	C	C	G	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	G	A	G	T	A	T	
	B-350101	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-350102	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-3507	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-3510	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-3511	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-3521	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-3524	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-3529	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-3540N	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-3541	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-3542	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-5305	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	A	T	
	B-3523	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	C	G	T	T	T	
	B-3546	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	T	C	T	T	C	A	G	A	G	G	A	G	G	A	T	G	G	A	T	G	T	A	T	
	B-5801	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	A	G	G	A	A	G	A	T	G	G	A	C	G	T	A	T	
	B-5804	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	A	G	G	A	A	G	A	T	G	G	A	C	G	T	A	T	
	B-5809	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	A	G	G	A	A	G	A	T	G	G	A	C	G	T	A	T	
	B-5802	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	T	G	G	A	A	G	A	T	G	G	A	C	G	T	A	T	
	B-5806	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	T	G	G	A	A	G	A	T	G	G	A	C	G	T	A	T	
	B-5807	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	T	G	G	A	A	G	A	T	G	G	A	C	G	T	A	T	
	B-570101	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	G	T	G	A	A	A	G	A	T	G	G	A	T	G	T	A	T
	B-5706	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	G	T	G	A	A	A	G	A	T	G	G	A	T	G	T	A	T
	B-5708	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	G	T	G	A	A	A	G	A	T	G	G	A	T	G	T	A	T
	B-1310	A	G	C	T	C	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	T	G	T	A	C		
	B-1540	A	G	C	T	C	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	T	G	T	A	C		
	B-1810	A	G	C	T	C	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	C	G	T	A	C		
	B-4003	A	G	C	T	C	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-4020	A	G	C	T	C	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-1552	A	G	C	T	C	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	G	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-4038	A	G	C	T	C	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-3515	A	G	C	T	C	C	G	C	C	G	A	G	T	G	T	A	T	G	T	C	T	T	C	A	G	A	G	G	A	A	G	A	T	G	G	A	C	G	T	A	T		
	B-570102	A	G	C	T	C	C	G	C	C	C	C	T	G	T	G	C	A	T	G	A	C	A	T	C	A	G	G	T	G	A	A	A	G	A	T	G	G	A	T	G	T	A	T
	B-1306	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	G	G	A	A	G	A	T	G	G	A	T	G	T	A	C		
	B-1542	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	A	G	A	T	G	G	A	T	G	T	A	C		
	B-5512	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-5501	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-5502	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-5505	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-5510	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-5516	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-5401	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	G	G	T	A	C		
	B-5402	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	G	G	T	A	C		
	B-5507	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	G	G	T	A	C		
	B-5503	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	T	G	A	G	A	T	G	G	A	A	G	T	A	C	
	B-3917	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-5610	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-5511	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	A	C	A	G	A	T	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-5901	A	G	C	T	T	A	G	C	C	A	C	G	T	G	T	A	T	G	T	C	T	T	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-1303	A	G	C	T	T	A	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	A	G	A	T	G	G	A	T	G	T	A	C		
	B-4410	A	G	C	T	T	A	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	T	T		
	B-4901	A	G	C	T	T	A	G	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C		
	B-4902	A	G	C	T	T	A	G	C	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	A	G	A	T	G	G	A	A	G	T	A	C	
	B-4903	A	G	C	T	T	A	G	C	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	A	G	A	T	G	G	A	C	G	T	A	C	
	B-5001	A	G	C	T	T	A	G	C	C	C	C	T	G	T	G	T	A	T	G	T	C	T	C	C	A	G	A	C	G	A	A												

5	B-5607	A	G	C	T	T	A	G	C	C	C	T	G	T	G	C	A	T	A	T	C	T	A	C	A	G	A	C	G	A	G	A	G	A	T	G	G	A	G	A	G	T	A	T	
	B-5602	A	G	C	T	T	A	G	C	C	C	C	T	G	T	G	C	A	T	A	T	C	T	A	C	A	G	A	C	G	A	G	A	G	A	T	G	G	A	G	A	G	T	A	C
	B-5604	A	G	C	T	T	A	G	C	C	C	C	T	G	T	G	C	A	T	A	T	C	T	A	C	A	G	A	C	G	A	G	A	G	A	T	G	G	A	G	A	G	T	A	C
	B-5608	A	G	C	T	T	A	G	C	C	C	C	T	G	T	G	C	A	T	A	A	G	T	A	C	A	G	A	C	G	A	G	A	G	A	T	G	G	A	G	A	G	T	A	T
	B-1302	A	G	C	T	T	A	G	C	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	G	A	T	G	G	G	A	T	G	T	A	T
	B-1308	A	G	C	T	T	A	G	C	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	G	A	T	G	G	G	A	T	G	T	A	T
	B-1309	A	G	C	T	T	A	G	C	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	G	A	T	G	G	G	A	T	G	T	A	T
	B-1301	A	G	C	T	T	A	G	C	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	G	A	T	G	G	G	A	T	G	T	A	T
	B-1307N	A	G	C	T	T	A	G	C	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	G	A	T	G	G	G	A	T	G	T	A	T
	B-1311	A	G	C	T	T	A	G	C	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	G	A	T	G	G	G	A	T	G	T	A	C
	B-4048	A	G	C	T	T	A	G	C	C	C	G	A	G	T	G	C	A	T	A	T	C	T	C	C	A	G	A	C	G	A	G	A	G	A	T	G	G	G	A	A	G	T	A	C
	B-1401	A	G	C	T	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	T	C	T	G	C	A	G	T	G	G	A	G	A	G	A	T	G	G	A	G	A	G	T	A	T
	B-1402	A	G	C	T	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	T	C	T	G	C	A	G	T	G	G	A	G	A	G	A	T	G	G	A	G	A	G	T	A	T
	B-1404	A	G	C	T	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	A	C	T	G	C	A	G	T	G	G	A	G	A	G	A	T	G	G	A	G	A	G	T	A	T
	B-1405	A	G	C	T	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	T	C	T	G	C	A	G	A	C	G	A	G	A	G	A	T	G	G	A	G	A	G	T	A	C
	B-140601	A	G	C	T	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	T	C	T	G	C	A	G	A	C	G	A	G	A	G	A	T	G	G	A	G	A	G	T	A	C
	B-140602	A	G	C	T	T	C	G	C	C	C	A	C	G	T	G	T	A	T	A	T	C	T	G	C	A	G	A	C	G	A	G	A	G	A	T	G	G	A	G	A	G	T	A	T
	B-3916	A	G																																										

B-5307	A G C T	T C G C C	C T G T	G T A T A	T C T T	C A G A G	G A G A A	G A T G	G G A C	G T C T
B-3503	A G C T	T C G C C	C T G T	G C A T G	T C T T	C A G A G	G A G A G	G A T G	G G A C	G T A T
B-3513	A G C T	T C G C C	C T G T	G C A T G	T C T T	C A G A G	G A G A G	G A T G	G G A C	G T A T
B-3536	A G C T	T C G C C	C T G T	G C A T G	T C T T	C A G A G	G A G A G	G A T G	G G A C	G T A T
B-5304	A G C T	T C G C C	C T G T	G C A T G	T C T T	C A G A G	G A G A A	G A T G	G G A C	G T A T
B-5611	A G C T	T C G C C	C T G T	G C A T G	T C T A	C A G A G	G A G A G	G A T G	G A G A	G T A T
B-3533	A G C T	T C G C C	G A G T	G C A T G	T C T T	C A G A G	G A G A G	G A T G	G G A C	G T A T
B-4036	A G C T	T C G C C	G A G T	G C A T A	T C T C	C A G A G	G A G A G	G A T G	G G A A	G T A C
B-4807	A G C T	T C G C C	G A G T	G C A T A	T C T C	C A G A G	G A G A G	C A T G	G A G A	G T A C
B-7301	A G C T	T C G C C	G A G T	G T A T A	T C T G	C A G A C	G T G G G	G A T G	G A G A	G T A T

10

15

20

25

30

TABLE XII

[illegible]

	DRB1-0429	ATGAGAGTACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-0430	ATGAGAGTACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-0445	ATGAGAGTACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-0448	ATGAGAGTACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
5	DRB1-0431	ATGAGAGTACGTAGTACGAGGACCAGAGTGGTGGGT	TGGT
	DRB1-0424	ATGAGAGTACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-0425	ATGAGAGTACGTAGTACGAGGACTACAGTGGTGGGT	TGTG
	DRB1-0436	ATGAGAGTACGTAGTACGAGGACTACAGCGGTGGGT	TGTG
10	DRB1-0447	ATGAGAGTACGTAGTACGAGGACTACAGCGGTGGGT	TGGT
	DRB1-0415	ATGAGAGTACGTAGTAGGAGGACTACAGCGGTGGGT	TGTG
	DRB1-040302	ATGAGAGTACGTAGTAGGAGGACCAGAGAGGTGGGT	TGTG
	DRB1-0435	ATGAGAGTACGTAGTTCGAGGACCAGAACGGTGGGT	TGGT
15	DRB1-0442	ATGAGAGTACGTAGTTCGAGGACCAGAGCGGTGGGT	TGTG
	DRB1-0428	ATGAGAGTACGTAGTTCGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-0443	ATGAGAGTACGTAGTTCGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-1122	ATGAGAGTACGTAGTTCGAGGACTACAGCGGTGGGT	TGGT
20	DRB1-0406	ATGAGAGTCCGTAGTACGAGGACCAGAGAGGTGGGT	TGTG
	DRB1-0446	ATGAGAGTCCGTAGTACGAGGACCAGAGAGGTGGGT	TGTG
	DRB1-0420	ATGAGAGTCCGTAGTACGAGGACCAGAGAGGTGGGT	TGGT
	DRB1-0421	ATGAGAGTCCGTAGTACGAGGACCAGAACGGTGGGT	TGGT
25	DRB1-0419	ATGAGAGTCCGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-1410	ATGAGAGTTCGTAGTAGGAGGACCAGAGAGGTGGGT	TGTG
	DRB1-1332	CTGAGAGAACGTAGTACGAGGACAAACGACGGTGGGT	TGTG
	DRB1-1340	CTGAGAGAACGTAGTACGAGGACAAACGACGGTGGGT	TGTG
30	DRB1-1353	CTGAGAGAACGTAGTACGAGGACAAACGACGGTGGGT	TGTG
	DRB1-1336	CTGAGAGAACGTAGTACGAGGACAAACGACGGTGGGT	TGGT
	DRB1-1424	CTGAGAGAACGTAGTACGAGGACAAAGGCCTGGGTGGGT	TGGT
	DRB1-030201	CTGAGAGAACGTAGTACGAGGACCAGAACGGGTGGGT	TGGT
	DRB1-030202	CTGAGAGAACGTAGTACGAGGACCAGAACGGGTGGGT	TGGT
	DRB1-0303	CTGAGAGAACGTAGTACGAGGACCAGAACGGGTGGGT	TGTG
	DRB1-0306	CTGAGAGAACGTAGTACGAGGACCAGAACGGGTGGGT	TGTG
	DRB1-1419	CTGAGAGAACGTAGTACGAGGACCAGAACGGGTGGGT	TGGT
	DRB1-1429	CTGAGAGAACGTAGTACGAGGACCAGAGCGGTGGGT	CTGTG
	DRB1-1406	CTGAGAGAACGTAGTACGAGGACCAGAGCGGTGGGT	TGTG
	DRB1-1402	CTGAGAGAACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-1409	CTGAGAGAACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-1413	CTGAGAGAACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-1446	CTGAGAGAACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-1447	CTGAGAGAACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT
	DRB1-1448	CTGAGAGAACGTAGTACGAGGACCAGAGCGGTGGGT	TGGT



5	DRB1-1403	CTGAGAGAAACGT	AGTACGAGGACCACAGTGGT	GGGTTGGT
	DRB1-140302	CTGAGAGAAACGT	AGTACGAGGACCACAGTGGT	GGGTTGGT
	DRB1-1412	CTGAGAGAAACGT	AGTACGAGGACCACAGTGGT	GGGTTGTG
	DRB1-1418	CTGAGAGAAACGT	AGTATGAGGACCGGAGAGGT	GGGTTGTG
	DRB1-1326	CTGAGAGAAACGT	AGTATGAGGACTACAGCGGT	GGGTTGGT
10	DRB1-1427	CTGAGAGAAACGT	AGTACGAGGACTACAGTGGT	GGGTTGGT
	DRB1-1334	CTGAGAGAAACCT	AGTTCGAGGACAACGACGGT	GGGTTGGT
	DRB1-0319	CTGAGAGAAACGT	AGTTCGAGGACAAGAAGGGT	GGGTTGTG
	DRB1-1310	CTGAGAGAAACGT	AGTTCGAGGACAACAACGGT	GGGTTGTG
	DRB1-130101	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGTG
	DRB1-130102	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGTG
	DRB1-130103	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGTG
	DRB1-1315	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGTG
	DRB1-1327	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGTG
	DRB1-1328	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGTG
15	DRB1-1335	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGTG
	DRB1-1351	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGTG
	DRB1-1359	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGTG
	DRB1-1361	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGTG
	DRB1-1316	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGAT
	DRB1-130201	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGGT
	DRB1-130202	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGGT
	DRB1-1331	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGGT
	DRB1-1339	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGGT
	DRB1-1341	CTGAGAGAAACGT	AGTTCGAGGACAACGACGGT	GGGTTGGT
20	DRB1-1309	CTGAGAGAAACGT	AGTTCGAGGACAAGGCCGGT	GGGTTGTG
	DRB1-1306	CTGAGAGAAACGT	AGTTCGAGGACAACAGCGGT	GGGTTGTG
	DRB1-1356	CTGAGAGAAACGT	AGTTCGAGGACCACAGCGGT	GGGTTGGT
	DRB1-0311	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-0324	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-0320	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGCTGTG
	DRB1-030101	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
25	DRB1-030102	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-0307	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-0312	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-0313	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-0315	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-0316	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-0318	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-0322	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-0323	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGTG
	DRB1-030501	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGGT
30	DRB1-030502	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGGT
	DRB1-0309	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGGT
	DRB1-0314	CTGAGAGAAACGT	AGTTCGAGGACCAGAAAGGT	GGGTTGGT
	DRB1-1421	CTGAGAGAAACGT	AGTTCGAGGACCAGAAACGGT	GGGTTGTG

5	DRB1-1417	CTGAGAGAAACGT	AGTT	CGAGGACC	CAGAGCGGT	GGGT	TGTG
	DRB1-1430	CTGAGAGAAACGT	AGTT	CGAGGACC	CAGAGCGGT	GGGT	TGGT
	DRB1-1433	CTGAGAGAAACGT	AGTT	CGAGGACC	CAGAGAGGT	GGGT	TGTG
	DRB1-1320	CTGAGAGAAACGT	AGTT	CGAGGACC	ACGACGGT	GGGT	TGTG
	DRB1-1329	CTGAGAGAAACGT	AGTT	CGAGGACC	ACGACGGT	GGGT	TGGT
10	DRB1-1342	CTGAGAGAAACGT	AGTT	CGAGGACT	TACAGCGGT	GGGT	TGTG
	DRB1-1305	CTGAGAGAAACGT	AGTT	CGAGGACT	TACAGCGGT	GGGT	TGGT
	DRB1-1350	CTGAGAGAAACGT	AGTT	CGAGGACT	TACAGCGGT	GGGT	TGGT
	DRB1-1318	CTGAGAGAAACGT	AGTT	CGAGGACT	TACAGTGGT	GGGT	TGTG
	DRB1-1116	CTGAGAGAAACGT	AGTT	GGAGGACA	ACGACGGT	GGGT	TGTG
15	DRB1-1120	CTGAGAGAAACGT	AGTT	GGAGGACA	ACGACGGT	GGGT	TGGT
	DRB1-0308	CTGAGAGAAACGT	AGTT	GGAGGACC	AGAAGGGT	GGGT	TGTG
	DRB1-0310	CTGAGAGAAACGT	AGTT	GGAGGACC	AGAAGGGT	GGGT	TGTG
	DRB1-1343	CTGAGAGAAACGT	AGTT	GGAGGACC	ACGACGGT	GGGT	TGTG
	DRB1-1109	CTGAGAGAAACGT	AGTT	GGAGGACT	TACAGCGGT	GGGT	TGGT
20	DRB1-1128	CTGAGAGAAACGT	AGTT	GGAGGACT	TACAGCGGT	GGGT	TGGT
	DRB1-1140	CTGAGAGAAACGT	AGTT	GGAGGACT	ACGACGGT	GGGT	TGTG
	DRB1-1115	CTGAGAGGACTT	AGTT	GGAGGACT	TACAGCGGT	GGGT	TGGT
	DRB1-1124	CTGAGAGGACGT	AGTT	GGAGGACT	TACAGCGGT	GGGT	TGGT
	DRB1-1362	CTGAGAGGACTT	AGTT	CGAGGACT	TACAGCGGT	GGGT	TGGT
25	DRB1-1144	CTGAGAGTACGC	AGTT	GGAGGACT	TACAGCGGT	GGGT	TGTG
	DRB1-130301	CTGAGAGTACGT	AGTAC	CGAGGACA	ACAACGGT	GGGT	TGGT
	DRB1-130302	CTGAGAGTACGT	AGTAC	CGAGGACA	ACAACGGT	GGGT	TGGT
	DRB1-1333	CTGAGAGTACGT	AGTAC	CGAGGACA	ACAACGGT	GGGT	TGGT
	DRB1-1337	CTGAGAGTACGT	AGTAC	CGAGGACA	ACAACGGT	GGGT	TGGT
30	DRB1-1338	CTGAGAGTACGT	AGTAC	CGAGGACA	ACGACGGT	GGGT	TGGT
	DRB1-1312	CTGAGAGTACGT	AGTAC	CGAGGACA	ACAGCGGT	GGGT	TGGT
	DRB1-1313	CTGAGAGTACGT	AGTAC	CGAGGACA	ACAGTGGT	GGGT	TGGT
	DRB1-1348	CTGAGAGTACGT	AGTAC	CGAGGACA	ACGACGGT	GGGT	TGTG
	DRB1-1358	CTGAGAGTACGT	AGTAC	CGAGGACA	ACAGCGGT	GGGT	CTGTG
	DRB1-0317	CTGAGAGTACGT	AGTAC	CGAGGACC	AGAAGGGT	GGGT	TGGT
	DRB1-0434	CTGAGAGTACGT	AGTAC	CGAGGACC	AGAACGGT	GGGT	TGGT
	DRB1-0820	CTGAGAGTACGT	AGTAC	CGAGGACT	TACAGTGGT	GGGT	TGTG
	DRB1-130701	CTGAGAGTACGT	AGTAC	CGAGGACT	TACAGCGGT	GGGT	TGGT
	DRB1-1349	CTGAGAGTACGT	AGTAC	CGAGGACT	TACAGCGGT	GGGT	TGGT
	DRB1-1347	CTGAGAGTACGT	AGTAC	CGAGGACT	TACAGTGGT	GGGT	TGGT
	DRB1-1355	CTGAGAGTACGT	AGTAC	CGAGGACT	TACAGTGGT	GGGT	TGGT

5	DRB1-1141	CTGAGAGTACGTAGTAGGAGGACTACGACGGTGGGTTGTG
	DRB1-1137	CTGAGAGTACGTAGTAGGAGGACTACAGCGGTGGGTTGGT
	DRB1-1425	CTGAGAGTACGTAGTAGGAGGACTACAGCGGTGGGTTGGT
	DRB1-130702	CTGAGAGTACGTAGTAGGAGGACTACAGCGGTGGGTTGGT
	DRB1-1442	CTGAGAGTACGTAGTTCGAGGACCGGAGAGGTGGGTTGGT
10	DRB1-1304	CTGAGAGTACGTAGTTCGAGGACAACGACGGTGGGTTGTG
	DRB1-1322	CTGAGAGTACGTAGTTCGAGGACAACGACGGTGGGTTGTG
	DRB1-1352	CTGAGAGTACGTAGTTCGAGGACAACGACGGTGGGTTGTG
	DRB1-1323	CTGAGAGTACGTAGTTCGAGGACAACGACGGTGGGTTGGT
	DRB1-1324	CTGAGAGTACGTAGTTCGAGGACTACGACGGTGGGTTGTG
15	DRB1-1354	CTGAGAGTACGTAGTTCGAGGACTACGACGGTGGGTTGTG
	DRB1-1311	CTGAGAGTACGTAGTTCGAGGACTACAGCGGTGGGTTGTG
	DRB1-1330	CTGAGAGTACGTAGTTCGAGGACAACAGCGGTGGGTTGGT
	DRB1-1325	CTGAGAGTACGTAGTTCGAGGACCACAGCGGTGGGTTGGT
	DRB1-131401	CTGAGAGTACGTAGTTCGAGGACTACAGCGGTGGGTTGGT
20	DRB1-1321	CTGAGAGTACGTAGTTCGAGGACTACAGCGGTGGGTTGGT
	DRB1-1346	CTGAGAGTACGTAGTTCGAGGACTACAGCGGTGGGTTGGT
	DRB1-1344	CTGAGAGTACGTAGTTCGAGGACCAGAGCGGTGGGTTGTG
	DRB1-0325	CTGAGAGTACGTAGTTCGAGGACCAGAAGGGTGGGTTGTG
	DRB1-1102	CTGAGAGTACGTAGTTCGAGGACAACGACGGTGGGTTGTG
25	DRB1-1121	CTGAGAGTACGTAGTTCGAGGACAACGACGGTGGGCTGTG
	DRB1-1118	CTGAGAGTACGTAGTTCGAGGACAACAGCGGTGGGTTGTG
	DRB1-1114	CTGAGAGTACGTAGTTCGAGGACAACGACGGTGGGTTGGT
	DRB1-1345	CTGAGAGTACGTAGTTCGAGGACAACGACGGTGGGTTGGT
	DRB1-1119	CTGAGAGTACGTAGTTCGAGGACAACAGCGGTGGGTTGGT
30	DRB1-1131	CTGAGAGTACGTAGTTCGAGGACAACAGCGGTGGGTTGGT
	DRB1-1145	CTGAGAGTACGTAGTTCGAGGACAACAGTGGTGGGTTGGT
	DRB1-1136	CTGAGAGTACGTAGTTCGAGGACCACGACGGTGGGTTGTG
	DRB1-1107	CTGAGAGTACGTAGTTCGAGGACCAGAAGGGTGGGTTGTG
	DRB1-1142	CTGAGAGTACGTAGTTCGAGGACCACAGCGGTGGGTTGTG
	DRB1-1134	CTGAGAGTACGTAGTTCGAGGACCAGAGCGGTGGGTTGTG
	DRB1-110801	CTGAGAGTACGTAGTTCGAGGACCACAGCGGTGGGTTGGT
	DRB1-110802	CTGAGAGTACGTAGTTCGAGGACCACAGCGGTGGGTTGGT
	DRB1-1126	CTGAGAGTACGTAGTTCGAGGACCAGAGCGGTGGGTTGGT
	DRB1-1103	CTGAGAGTACGTAGTTCGAGGACTACGACGGTGGGTTGTG
	DRB1-110601	CTGAGAGTACGTAGTTCGAGGACTACAGCGGTGGGCTGTG
	DRB1-110602	CTGAGAGTACGTAGTTCGAGGACTACAGCGGTGGGCTGTG

	DRB1-1135	CTGAGAGTACGTAGTTGGACGACTACAGCGGTGGGTGTGTG
	DRB1-110401	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGTG
	DRB1-110402	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGTG
	DRB1-1143	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGTG
	DRB1-1146	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGTG
	DRB1-1138	CTGAGAGTACGTAGTTGGGGGACTACAGCGGTGGGTGTGTG
5	DRB1-1125	CTGAGAGTACGTAGTTGGAGGACTACAGTGGTGGGTGTGTG
	DRB1-1111	CTGAGAGTACGTAGTTGGAGGACTACGACGGTGGGTGTGGT
	DRB1-1133	CTGAGAGTACGTAGTTGGACGACTACAGCGGTGGGTGTGGT
	DRB1-110101	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGGT
	DRB1-110102	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGGT
	DRB1-110103	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGGT
10	DRB1-110104	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGGT
	DRB1-110105	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGGT
	DRB1-112701	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGGT
	DRB1-112702	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGGT
	DRB1-1130	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGGT
	DRB1-1139	CTGAGAGTACGTAGTTGGAGGACTACAGCGGTGGGTGTGGT
	DRB1-1123	CTGAGAGTACGTAGTTGGAGGACTACAGTGGTGGGTGTGGT
	DRB1-1132	CTGAGAGTACGTAGTTGGAGGACTACAGTGGTGGGTGTGGT
15	DRB1-131402	CTGAGAGTACGTAGTTTGAGGACTACAGCGGTGGGTGTGGT
	DRB1-0304	CTGAGAGTCCGTAGTTGAGGACCAGAAGGGTGGGTGTGTG
	DRB1-1129	CTGAGAGTCCGTAGTTGGAGGACTACAGCGGTGGGTGTGGT
	DRB1-1147	CTGAGAGTCCGTAGTTGGAGGACTACAGCGGTGGGTGTGTG
	DRB1-1360	CTGAGAGTCCGTAGTATGAGGACCACAGCGGTGGGTGTGGT
20	DRB1-1441	CTGAGAGTTCCGTAGTACGAGGACCAGAGCGGTGGGTGTGGT
	DRB1-1308	CTGAGAGTTCCGTAGTACGAGGACAACGACGGTGGGTGTGTG
	DRB1-1319	CTGAGAGTTCCGTAGTACGAGGACAACGACGGTGGGTGTGTG
	DRB1-140502	CTGAGAGTTCCGTAGTACGAGGACC GGAGAGGTGGGTGTGTG
	DRB1-1423	CTGAGAGTTCCGTAGTACGAGGACC GGAGAGGTGGGTGTGTG
	DRB1-1420	CTGAGAGTTCCGTAGTACGAGGACCAGAGCGGTGGGTGTGTG
	DRB1-1357	CTGAGAGTTCCGTAGTTGAGGACAACGACGGTGGGTGTGTG
25	DRB1-0321	CTGAGAGTTCCGTAGTTGAGGACCAGAAGGGTGGGTGTGTG
	DRB1-1416	CTGAGAGTTCCGTAGTAGGAGGACAACGACGGTGGGTGTGTG
	DRB1-1117	CTGAGAGTTCCGTAGTAGGAGGACC GGAGAGGTGGGTGTGTG
	DRB1-140101	CTGAGAGTTCCGTAGTAGGAGGACC GGAGAGGTGGGTGTGTG
	DRB1-140102	CTGAGAGTTCCGTAGTAGGAGGACC GGAGAGGTGGGTGTGTG
	DRB1-1408	CTGAGAGTTCCGTAGTAGGAGGACC GGAGAGGTGGGTGTGTG
	DRB1-1426	CTGAGAGTTCCGTAGTAGGAGGACC GGAGAGGTGGGTGTGTG
30	DRB1-1438	CTGAGAGTTCCGTAGTAGGAGGACC GGAGAGGTGGGTGTGTG
	DRB1-1439	CTGAGAGTTCCGTAGTAGGAGGACC GGAGAGGTGGGTGTGTG
	DRB1-1432	CTGAGAGTTCCGTAGTAGGAGGACC GGAGCGGTGGGTGTGTG
	DRB1-1434	CTGAGAGTTCCGTAGTAGGAGGACC GGAGCGGTGGGTGTGTG

5	DRB1-1113	CTGAGAGTTCGTAGTTGGAGGACC	GGAGCGGTGGGTTGTG
	DRB1-1435	CTGAGAGTTCGTAGTTGGAGGACC	GGAGAGGTGGGTTGTG
	DRB1-1437	CTGAGAGTTCGTAGTATGAGGACA	AGGCCGGTGGGTTGTG
	DRB1-1445	CTGAGAGTTCGTAGTATGAGGACA	GGAGAGGTGGGTTGTG
	DRB1-140501 DRB1-1443	CTGAGAGTTCGTAGTATGAGGACC	GGAGAGGTGGGTTGTG
10	DRB1-1110 DRB1-111201 DRB1-111202	CTGAGAGTTCGTAGTTGGAGGACT	ACAGCGGTGGGTTGGT
	DRB1-1414 DRB1-1436	CTGAGAGTTCGTAGTACGAGGACC	GGAGAGGTGGGTTGGT
	DRB1-140701 DRB1-140702	CTGAGAGTTCGTAGTAGGAGGACC	GGAGAGGTGGGTTGGT
	DRB1-1422	CTGAGAGTTCGTAGTAGGAGGACT	ACAGCGGTGGGTTGGT
	DRB1-1440 DRB1-1444	CTGAGAGTTCGTAGTACGAGGACC	GGAGAGGTGGGTTGGT
15	DRB1-120101 DRB1-120102 DRB1-1206 DRB1-1207 DRB1-1208 DRB1-1209	GTGAGAGCTCCTAGTTTGAGGACA	ACAGCGGTGGGCTGTG
	DRB1-120302	GTGAGAGCTCCTAGTTTGAGGACA	ACAGCGGTGGGTTGTG
	DRB1-1204	GTGAGAGCTCCTAGTTTGAGGACA	ACAGCGGTGGGCTGTG
	DRB1-120201 DRB1-120202	GTGAGAGCTCCTAGTTTGAGGACT	ACAGCGGTGGGCTGTG
	DRB1-0816 DRB1-0818 DRB1-0825	GTGAGAGTACGTAGTACGAGGACT	ACAGTGGTGGGTTGGT
20	DRB1-0810 DRB1-0812	GTGAGAGTACGTAGTACGAGGACA	ACAGTGGTGGGTTGTG
	DRB1-080302 DRB1-0814 DRB1-0819 DRB1-0823	GTGAGAGTACGTAGTACGAGGACA	ACAGTGGTGGGTTGGT
	DRB1-0813	GTGAGAGTACGTAGTACGAGGACC	ACAGTGGTGGGTTGGT
	DRB1-080401 DRB1-080404 DRB1-0806	GTGAGAGTACGTAGTACGAGGACT	ACAGTGGTGGGTTGTG
	DRB1-0822	GTGAGAGTACGTAGTACGAGGACT	ACAGTGGTGGGCTGTG
30	DRB1-0805 DRB1-0824	GTGAGAGTACGTAGTACGAGGACT	ACAGCGGTGGGTTGGT

5	DRB1-080101	GTGAGAGTACGTAGTACGAGGACTACAGTGGTGGGT	TTGGT
	DRB1-080102	GTGAGAGTACGTAGTACGAGGACTACAGTGGTGGGT	TTGGT
	DRB1-080201	GTGAGAGTACGTAGTACGAGGACTACAGTGGTGGGT	TTGGT
	DRB1-080202	GTGAGAGTACGTAGTACGAGGACTACAGTGGTGGGT	TTGGT
	DRB1-080203	GTGAGAGTACGTAGTACGAGGACTACAGTGGTGGGT	TTGGT
	DRB1-0807	GTGAGAGTACGTAGTACGAGGACTACAGTGGTGGGT	TTGGT
	DRB1-0811	GTGAGAGTACGTAGTACGAGGACTACAGTGGTGGGT	TTGGT
10	DRB1-080402	GTGAGAGTACGTAGTACGAGGACTACAGTGGTGGGT	TTGTT
	DRB1-080403	GTGAGAGTACGTAGTACGAGGACTACAGTGGTGGGT	TTGTT
	DRB1-0808	GTGAGAGTACGTAGTAGGAGGACTACAGTGGTGGGT	TTGGT
	DRB1-0815	GTGAGAGTACGTAGTAGGAGGACAACAGTGGTGGGT	TTGGT
	DRB1-0817	GTGAGAGTACGTAGTTGAGGACTACAGTGGTGGGT	TTGGT
	DRB1-1317	GTGAGAGTACGTAGTTGAGGACAACGACGGTGGGT	TTGTG
	DRB1-1105	GTGAGAGTACGTAGTTGAGGACTACAGC	GGTGGGT
15	DRB1-0809	GTGAGAGTTCGTAGTACGAGGACTACAGTGGTGGGT	TTGGT
	DRB1-0821	GTGAGAGTTCGTAGTACGAGGACTACAGTGGTGGGT	TTGGT
	DRB1-1415	GTGAGAGTTCGTAGTACGAGGACTACAGTGGTGGGT	TTGTG
	DRB1-1205	GTGAGAGTTCCTAGTTCGAGGACAACAGC	GGTGGGT
	DRB1-1404	GTGAGAGTTCGTAGTAGGAGGACC	GGAGAGGTGGGT
	DRB1-1411	GTGAGAGTTCGTAGTAGGAGGACC	GGAGAGGTGGGT
	DRB1-1428	GTGAGAGTTCGTAGTAGGAGGACC	GGAGAGGTGGGT
20	DRB1-1431	GTGAGAGTTCGTAGTAGGAGGACC	GGAGC
	DRB1-1507	GGGAGAGTCCGTAGTATGAGGACAAGGCC	GGTGGGT
	DRB1-1511	GGGAGAGTCCGTAGTATGAGGACAAGGCC	GGTGGGT
	DRB1-1605	GGGAGAGTCCGTAGTATGAGGACAACAGC	GGTGGGT
	DRB1-1607	GGGAGAGTCCGTAGTATGAGGACAACAGC	GGTGGGT
	DRB1-160201	GGGAGAGTCCGTAGTATGAGGACC	ACAGC
	DRB1-160202	GGGAGAGTCCGTAGTATGAGGACC	ACAGC
25	DRB1-160101	GGGAGAGTCCGTAGTATGAGGACTACAGC	GGTGGGT
	DRB1-160102	GGGAGAGTCCGTAGTATGAGGACTACAGC	GGTGGGT
	DRB1-1603	GGGAGAGTCCGTAGTATGAGGACTACAGC	GGTGGGT
	DRB1-1604	GGGAGAGTCCGTAGTATGAGGACTACAGT	GGTGGGT
	DRB1-150104	GGGAGAGTCCGTAGTTTCGAGGACAAGGCC	GGTGGGT
	DRB1-1512	GGGAGAGTCCGTAGTTTCGAGGACAAGGCC	GGTGGGT
	DRB1-150202	GGGAGAGTCCGTAGTTTCGAGGACAAGGCC	GGTGGGT
30	DRB1-1510	GGGAGAGTCCGTAGTTTCGAGGACAACGAC	GGTGGGT
	DRB1-1508	GGGAGAGTCCGTAGTTTCGAGGACAAGGCC	GGTGGGT
	DRB1-150102	GGGAGAGTCCGTAGTTTCGAGGACAAGGCC	GGTGGGT
	DRB1-150101	GGGAGAGTCCGTAGTTTCGAGGACAAGGCC	GGTGGGT
	DRB1-150103	GGGAGAGTCCGTAGTTTCGAGGACAAGGCC	GGTGGGT

5	DRB1-150105	G	G	G	A	G	A	G	T	C	C	G	T	A	G	T	T	T	G	A	G	G	A	C	A	A	G	G	C	C	G	G	T	G	G	G	T	T	G	T	G		
	DRB1-1503	G	G	G	A	G	A	G	T	C	C	G	T	A	G	T	T	T	G	A	G	G	A	C	A	A	G	G	C	C	G	G	T	G	G	G	T	T	G	T	G		
	DRB1-1506	G	G	G	A	G	A	G	T	C	C	G	T	A	G	T	T	T	G	A	G	G	A	C	A	A	G	G	C	C	G	G	T	G	G	G	T	T	G	T	G		
	DRB1-1509	G	G	G	A	G	A	G	T	C	C	G	T	A	G	T	T	T	G	A	G	G	A	C	A	A	G	G	C	C	G	G	T	G	G	G	T	T	G	T	G		
	DRB1-1513	G	G	G	A	G	A	G	T	C	C	G	T	A	G	T	T	T	G	A	G	G	A	C	A	A	G	G	C	C	G	G	T	G	G	G	T	T	G	T	G		
10	DRB1-150201	G	G	G	A	G	A	G	T	C	C	G	T	A	G	T	T	T	G	A	G	G	A	C	A	A	G	G	C	C	G	G	T	G	G	G	T	T	G	G	T		
	DRB1-150203	G	G	G	A	G	A	G	T	C	C	G	T	A	G	T	T	T	G	A	G	G	A	C	A	A	G	G	C	C	G	G	T	G	G	G	T	T	G	G	T		
	DRB1-1505	G	G	G	A	G	A	G	T	C	C	G	T	A	G	T	T	T	G	A	G	G	A	C	A	A	G	G	C	C	G	G	T	G	G	G	T	T	G	T	G		
	DRB1-1504	G	G	G	A	G	A	G	T	C	C	G	T	A	G	T	T	T	G	A	G	G	A	C	T	A	A	G	G	C	C	G	G	T	G	G	G	T	T	G	T	G	
	DRB1-1608	G	G	G	A	G	A	G	A	A	C	G	T	A	G	T	A	T	G	A	G	G	A	C	T	A	C	A	G	C	C	G	G	T	G	G	G	T	T	G	G	T	
15	DRB1-090102	T	T	G	A	G	A	G	A	A	C	G	T	A	G	T	A	C	G	A	G	G	A	C	T	G	G	A	G	A	G	G	T	G	G	G	T	T	G	G	T		
	DRB1-0902	T	T	G	A	G	A	G	A	A	C	G	T	A	G	T	A	T	G	A	G	G	A	C	T	G	G	A	G	A	G	G	T	G	G	G	T	T	G	G	T		
	DRB1-010102	T	T	G	A	G	A	G	A	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	A	G	C	C	G	G	T	G	G	G	T	T	G	G	T
	DRB1-0108	T	T	G	A	G	A	G	T	A	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	A	G	C	C	G	G	T	G	G	G	T	T	G	G	T	
	DRB1-100101	T	T	G	A	G	A	G	T	A	C	G	C	A	G	T	A	C	G	A	G	G	A	C	C	G	G	A	G	C	C	G	G	T	G	G	G	T	T	G	G	T	
20	DRB1-100102	T	T	G	A	G	A	G	T	A	C	G	C	A	G	T	A	C	G	A	G	G	A	C	C	G	G	A	G	C	C	G	G	T	G	G	G	T	T	G	G	T	
	DRB1-0103	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	A	A	C	G	A	C	C	G	G	T	G	G	G	T	T	G	G	T	
	DRB1-0110	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	A	C	C	G	G	T	G	G	G	T	T	G	G	T		
	DRB1-0106	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	G	C	C	G	G	T	G	G	G	T	T	G	T	G		
	DRB1-0109	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	G	C	C	G	G	T	G	G	G	T	T	G	G	T		
25	DRB1-010202	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	A	G	C	C	G	T	G	G	G	C	T	G	T	G		
	DRB1-010201	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	A	G	C	C	G	T	G	G	G	C	T	G	T	G		
	DRB1-0104	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	A	G	C	C	G	T	G	G	G	T	T	G	T	G		
	DRB1-010101	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	A	G	C	C	G	G	T	G	G	G	T	T	G	G	T	
	DRB1-0105	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	A	G	C	C	G	G	T	G	G	G	T	T	G	G	T	
30	DRB1-0107	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	A	G	C	C	G	G	T	G	G	G	T	T	G	G	T	
	DRB1-0111	T	T	G	A	G	A	G	T	C	C	G	T	A	G	T	A	C	G	A	G	G	A	C	C	A	G	A	G	C	C	G	G	T	G	G	G	T	T	G	G	T	

General strategy for medium resolution typing is described below:

For medium resolution typing a maximally informative set of marker positions were determined. These consist of positions 98, 414, 539, 282, 571, 368, 256, 292, 238, 270, 453, 527, 502, 81, 268, 559, 92, 123 and 396 of HLA-A (numbering starts at the transcription start position of exon 1), positions 539, 419, 559, 412, 272, 362, 302, 363, 206, 369, 259, 97, 583, 292, 222, 527, 418, 435 and 571 of HLA-B (numbering starts at the transcription start position of exon 1), and positions 125, 196, 197, 227, 261, 286, 299, 308, 341 and 345 of HLA-DRB1 (numbering starts at the transcription start position of exon 1).

In general, the order of the positions is from the most informative to the least informative with respect to the selection criteria of frequent and rare HLA alleles (see list of frequent HLA alleles above). Thus the ten markers (HLA-A and HLA-B) that were selected for the fine typing strategy constitute the first ten markers of the set of 19 markers for the single pass classification into frequent and rare HLA alleles (HLA-A and HLA-B). Like with sequence-based HLA typing there are heterozygous combinations of HLA alleles that can not be resolved. However, there are fewer ambiguities with this method due to the mini-haplotypes that are provided.

Another object of the present invention is the use of said methodology of the invention is for screening of tissue donors, for example, bone marrow donors in registries for frequent and rare HLA types.

The description of the HLA alleles is based on the Anthony Nolan database ([www.ebi.ac.uk/imgt/hla/](http://www.ebi.ac.uk/imgt/hla/)).

In addition to the aforementioned method, the invention includes yet other arrangements which will emerge from the description that follows, which refers to examples of supports according to the invention, as well as the annexed figures and tables, wherein:



Figure 1 describes 19 positions covered by mini-haplotyping assays for discrimination of HLA-A mapped onto the HLA-A allele A\*010101 as reference. Black boxes indicate an extension position while grey boxes indicate polymorphisms that are captured by the annealing of the respective primer of the primer pool. Pools are used in forward and reverse. Numbering is according to the transcription start of the cDNA.

Figure 2 describes 19 positions covered by mini-haplotyping assays for discrimination of HLA-B mapped onto the HLA-B allele B\*070201 as reference. Black boxes indicate an extension position while grey boxes indicate polymorphisms that are captured by the annealing of the respective primer of the primer pool. Pools are used in forward and reverse. Numbering is according to the transcription start of the cDNA.

Figure 3 describes 10 positions covered by mini-haplotyping assays for discrimination of HLA-DRB1 mapped onto the HLA-DRB1 allele DRB1\*0101 as reference. Black boxes indicate an extension position while grey boxes indicate polymorphisms that are captured by the annealing of the respective primer of the primer pool. Pools are used in forward and reverse. Numbering is according to the transcription start of the cDNA.

Figure 4 describes 10 positions covered by mini-haplotyping assays for discrimination of HLA-A mapped onto the HLA-A allele A\*010101 as reference for the distinction of subgroups that can then be further analysed. Black boxes indicate an extension position while grey boxes indicate polymorphisms that are captured by the annealing of the respective primer of the primer pool. Pools are used in forward and reverse. Numbering is according to the transcription start of the cDNA.

Figure 5 describes 10 positions covered by mini-haplotyping assays for discrimination of HLA-B mapped onto the HLA-B allele B\*070201 as reference for the distinction of subgroups that can then be further analysed. Black boxes indicate

an extension position while grey boxes indicate polymorphisms that are captured by the annealing of the respective primer of the primer pool. Pools are used in forward and reverse. Numbering is according to the transcription start of the cDNA.

- 5 Figure 6 shows genotyping results of a CEPH family (1418, 01 = father, 02 = mother, 03 = child, 04 = child) for position HLA-B\_272. 1407,3 Da corresponds to the addition of C to primer 6, 7, 8, or 9; 1422,3 Da corresponds to the addition of T to primer 6, 7, 8, or 9; 1431,4 Da/ 1430,9 Da corresponds to the addition of A to primer 6, 7, 8, or 9; and 1447,4 Da/ 1448,5 Da corresponds to the addition of G to  
10 primer 6, 7, 8, or 9.

Table I represents HLA-A alleles captured by the 10 markers in the different subgroups and additional positions that have to be typed to resolve the subgroups.

- 15 Table II represents HLA-B alleles captured by the 10 markers in the different subgroups and additional positions that have to be typed to resolve the subgroups.

Table III represents HLA-DRB1 alleles captured by the 10 markers in the different subgroups and additional positions that have to be typed to resolve the subgroups.

20

- Table IV represents the list of the individual primers that are required to constitute the pools for mini-haplotyping of HLA-A (19 markers). The 10 markers required for the creation of subgroups are also contained. ^ refers to the base used to attach the mass/charge tag, CT refers to the mass difference of the mass/charge tag, sp  
25 means phosphorothioate group. The product analysed by mass spectrometry includes the base 5' of the most 5' phosphorothioate (sp).

- Table V represents the list of the individual primers that are required to constitute the pools for mini-haplotyping of HLA-B (19 markers). The 10 markers required  
30 for the creation of subgroups are also contained. ^ refers to the base used to attach the mass/charge tag, CT refers to the mass difference of the mass/charge tag, sp

means phosphorothioate group. The product analysed by mass spectrometry includes the base 5' of the most 5' sp.

5 Table VI represents the list of the individual primers that are required to constitute the pools for mini-haplotyping of HLA-DRB1 (10 markers). ^ refers to the base used to attach the mass/charge tag, CT refers to the mass difference of the mass/charge tag, sp means phosphorothioate group. The product analysed by mass spectrometry includes the base 5' of the most 5' sp.

10 Table VII represents the resolution that can be generated with the 19 markers for the distinction of the frequent HLA alleles in HLA-A.

Table VIII represents the resolution that can be generated with the 19 markers for the distinction of the frequent HLA alleles in HLA-B.

15

Table IX represents the resolution that can be generated with the 10 markers for the distinction of the frequent HLA alleles in HLA-DRB1.

20 Table X represents the list of HLA-A alleles that are resolved with the 10 markers for the creation of subgroups. Each subgroup is separated by an empty line. Frequent alleles are shaded in darker grey, while lighter grey indicates the position that primers are extended onto.

25 Table XI represents the list of HLA-B alleles that are resolved with the 10 markers for the creation of subgroups. Each subgroup is separated by an empty line. Frequent alleles are shaded in darker grey, while lighter grey indicates the position that primers are extended onto.

30 Table XII represents the list of HLA-DRB1 alleles that are resolved with the 10 markers for the creation of subgroups. Each subgroup is separated by an empty line. Frequent alleles are shaded in darker grey, while lighter grey indicates the position that primers are extended onto.

### Examples

#### 5 Example 1: Mini-haplotyping at position 272 of HLA-B by the modified GOOD-Assay

A locus specific PCR product of exon 2 and exon 3 of HLA-B is amplified with a set of primers published by the International Histocompatibility Working Group, Technical Manuals (Hurly, Fernandes-Vina, Gao, Middleton, Noreen, Ren and  
10 Smith; [www.ihwg.org/tmanual/Tmcontents.htm](http://www.ihwg.org/tmanual/Tmcontents.htm)). The PCR product is incubated with SAP to remove all excess dNTPs. Then a single base primer extension at position 272 in the PCR amplicon is carried out. The set of primers, to generate the mini-haplotypes is shown in Table V. Thereafter a 5'phosphodiesterase digest is applied to reduce the primers to a core sequence. After alkylation of the DNA  
15 backbone of the mini-haplotype fragments the products are transferred onto a MALDI target pre-coated with matrix. Alternatively the matrix solution can be mixed with the samples and transferred onto the MALDI target to dry. The MALDI target is introduced into a MALDI mass spectrometer and analysed. The mass spectra show one or two mass peaks and that correspond to specific mini-  
20 haplotypes.

#### PCR:

Forward primer, BAmp1 5'-G GGT CCC AGT TCT AAA GTC CCC ACG-3' (1.875 pmol), reverse primer, BAmp2 5'-CC ATC CCC GGC GAC CTA TAG  
25 GAG ATG-3' (1.875 pmol) and BAmp3 5'-AGG CCA TCC CGG CGG GCG ATC TAT-3' (1.875 pmol), 0.25 µl 10x PCR buffer (HiFi Platinum Taq), 0.3 µl MgSO<sub>4</sub> (50 mM), 0.2 µl of a mix of each dCTP, dATP, dGTP and dTTP (2 mM each), 0.25U engineered DNA polymerase (HiFi Platinum DNA Polymerase; Invitrogen) and 5 ng DNA fill to 3 µl with water. Cycling: 1. 94°C 3 min, 2. 94°C 20 sec, 3.  
30 64°C 30 sec, 4. 72°C 30 sec, steps 2 to 4 are repeated 35 times, 5. 72°C 5 min.

#### SAP digest:

1.75  $\mu$ l of 50 mM Tris-HCl and 0.25  $\mu$ l SAP (USB corporation, Cleveland, USA) are to add to the PCR product and this has to be incubated for 60 min at 37°C, followed by an incubation at 90°C for 10 min to denature the SAP enzyme.

5 Single Base Primer Extension:

To the SAP treated PCR product 2  $\mu$ l of an extension mix is to add. This mix contains 15 mM MgCl<sub>2</sub>, 0.1 mM of each of the four  $\alpha$ -S-ddNTPs, 5 pmol of the extension primers set and 0,4 U of Thermosequenase. Cycling: 1. 94°C 2 min, 2. 94°C 15 sec, 3. 58°C 20 sec, 4. 72°C 20 sec, steps 2 to 4 are repeated 50 times.

10

PDE digest:

To the extension product has to be added 0.5  $\mu$ l 0.5 M acetic acid and 1.5  $\mu$ l PDE (5.1U) and incubate for at lease 120 min at 37 °C.

15 Alkylation:

The alkylation is carried out by adding 21  $\mu$ l of an alkylation mix and incubate for 15 min at 40°C. Th alkylation mix contains 377 parts water free acetonitrile, 15 parts of 2M triethylamine/CO<sub>2</sub> (pH ~7.5), 75 parts 2mM Tris-HCl and 174 parts of methyl iodine.

20 The alkylation is to stopped by adding 10  $\mu$ l deionised water. 5  $\mu$ l of the resulting upper phase are to dilute in 10  $\mu$ l 40% acetonitrile.

For MALDI target preparation and measurement with the MALDI mass spectrometer 0.5  $\mu$ l of the final dilution are transferred onto a MALDI target pre-coated with matrix ( $\alpha$ -cyano-4-hydroxycinnamic acid methyl ester). Measurement  
25 was carried out in a Bruker Autoflex with typically -18 kV acceleration voltage, pulsed ion extraction with a delay of 200 ns, and detection in linear detection mode. Results for CEPH family 1418 are shown in figure 6.

30

Example 2: HLA-DR typing by the GOOD-Assay

A locus specific PCR for HLA-DRB is carried out. Therefore a set of allele-specific primers as listed below is used. These primers are published by J. Wu et al. in <http://www.ihwg.org/tmanual/TMcontents.htm> Chapter 10-B.

Name	Sequence
Amp1_DRB1_f20	5'-TTCTTGTGGSAGCTTAAGTT-3'
Amp2_DRB1_f21	5'-TTCCTGTGGCAGCCTAAGAGG-3'
Amp3_DRB1_f22	5'-CACGTTTCTTGGAGTACTCTAB-3'
Amp3-2_DRB1_f23	5'-CGTTTCTTGGAGTACTCTACGGG-3'
Amp3-3_DRB1_f23	5'-CGTTTCTTGGAGTACTCTACGTC-3'
Amp4_DRB1_f21	5'-GTTTCTTGGAGCAGGTTAAAC-3'
DR7_DRB1_f20	5'-CCTGTGGCAGGGTAARTATA-3'
DR9_DRB1_f18	5'-CCCAACCACGTTTCTTGA-3'
DR10_DRB1_f19	5'-AGACCACGTTTCTTGGAGG-3'
AmpB_DRB1_r18	5'-TCGCCGCTGCACYGTGAA-3'

5

This set of primers carries a high risk of co-amplifying genes for the other HLA-DRB chains, which results in unclear results. However, this is currently the best available option for the PCR of HLA-DRB1. In order to resolve the problem an additional mini-haplotyping test can be added. The mini-haplotyping assay HLA-DRB\_122-126 gives good resolution of HLA-DRB genes and allows the verification of results produced for typing of HLA-DRB1 PCR products. The identification of HLA-DRB1 genes is possible, as well as the identification of other amplified HLA-DRB genes which are present is possible. The set of primers listed below is used for the primer extension reaction. The details of the protocol are

15

identical to example 1.

Name	Sequence	CT	Masses				
			Primer	A	C	G	T
HLADR_1221_2f20	TGAAGAAATGACACTCAspTspG*spT	0	1487,5	-	-	-	1805,7
HLADR_1222_2f20	TGCAGAAATAGCACTCGspTspG*spT	0	1503,5	-	-	-	1821,7
HLADR_1223_2f20	TGAAGAAATGACACTCAspGspG*spT	0	1512,5	-	-	-	1830,7
HLADR_1224_2f20	TGAAGAAATGACACTTAspTspA*spT	0	1471,5	-	-	-	1789,7
HLADR_1225_2f20	TGAAGAAATGACACTCCspCspT*spC	-14	1510,6	-	-	-	1814,8
HLADR_1226_2f20	TGAAGAAATRCACTCAspCspC*spC	-28	1418,4	1717,7	1693,6	1733,7	-
HLADR_1227_2f20	TGAAGAAATGACACTCAspTspA*spC	-14	1456,5	-	-	-	1760,7
HLADR_1228_2f20	TGAAGAAWTGACACTCAspGspA*spC	0	1481,5	-	-	-	1799,7
HLADR_1229_2f20	TGAGGAAATGACACTCAspCspA*spC	-14	1441,5	-	-	1770,8	1745,7
HLADR_12210_2f20	TGAAGATATGACACTCAspCspA*spC	-14	1441,5	-	-	1770,8	1745,7
HLADR_12211_2f20	TGAAGAAATGACAYTCAspAspA*spC	0	1465,5	-	-	-	1783,7

Of the thirteen possible mini-haplotypes, four represent genes other than HLA-DRB1. The mini-haplotype GTGTT (1821.7 Da), AACAC in sense direction, represents with 100% certainty co-amplification of the HLA-DRB9 gene. The mini-haplotype ATACT (1760.8 Da), AGTAT in sense direction, represent either all  
5 HLA-DRB1\*07 alleles (except HLA-DRB1\*070102) or co-amplification of the HLA-DRB5 gene. The type TGTGT (1745.7 Da), AGTGT in sense direction, correspond to co-amplification or all variations of the HLA-DRB4 or HLA-DRB6 genes. Finally the type AGACT (1799.7 Da), AGTCT in sense direction, represent  
10 variants of HLA-DRB3 and HLA-DRB7 genes.

**This Page is Inserted by IFW Indexing and Scanning  
Operations and is not part of the Official Record**

**BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☒ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☒ **FADED TEXT OR DRAWING**
- ☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☒ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** \_\_\_\_\_

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.**